



SS1033-250

## May-Grunwald Giemsa Stain Kit

### Solutions Provided:

SSC1207-250 - **Jenner Stain Stock Solution**  
SSC1149-250 - **Giemsa Stain Stock Solution**  
SSC1208-250 - **Giemsa Buffer, pH 6.8**  
SSC1031-250 - **Acetic Acid, 1%**

### Working solution preparation (*prepare just before use*):

<b>Jenner Stain Working Solution</b>	<b>Giemsa Stain Working Solution</b>
Jenner Stain Stock Solution: 25mL	Giemsa Stain Stock Solution: 2.5mL
Giemsa Buffer, pH 6.8*: 25mL	Giemsa Buffer, pH 6.8*: 50mL

\*Final pH of working stain solutions will produce varying results. User preference may require the use of distilled/deionized water or the Giemsa Buffer, pH 7.2.

### Conventional Procedure

1. Deparaffinize and hydrate sections through alcohol to distilled water.
2. Place sections in 100% methanol for 3 minutes.
3. Place sections in a second change of 100% methanol for 3 minutes.
4. Stain sections in **Jenner Stain Working Solution** for **5-6 minutes**.
5. Stain sections in **Giemsa Stain Working Solution** for **45 minutes**.
6. Rinse slides briefly in distilled water.
7. Dip slides in **Acetic Acid, 1%** and check under microscope until desired differentiation has been achieved.
8. Rinse slides in distilled water
9. Dehydrate rapidly in 3 changes of 100% alcohol.
10. Clear and mount with appropriate mounting medium.

### Results:

Nuclei:	Blue
Cytoplasm:	Shades of pink
Bacteria:	Blue

### Notes:

1. After initial deparaffination and hydration: If sections have been fixed in Zenker, B-5 Solution or any other solution containing mercury be sure to prepare slides by placing in Gram's or Lugol's Iodine for 10 minutes, washing in tap water followed by Sodium Thiosulfate, 5% for 5 minutes. After one final wash in tap water and rinse in distilled water, proceed to **Step 2**.

