

**DANGER:**

Causes severe skin burns and serious eye damage. May cause damage to organs (eyes, skin, respiratory system) through prolonged or repeated exposure.



Deionized Water  
Formic Acid

7732-18-5  
64-18-6



*For In Vitro Diagnostic Use Only*

# Formic Acid, Decalcifier

# CalFor\*™

**DIRECTIONS FOR USE:**

Calfor\*™ is a fast, yet exceptionally gentle bone decalcifier that has been formulated from semiconductor grade formic acid and highly purified reagent grade water. Calfor\* is formulated as a general purpose decalcifier for all routine specimens, and because it will safely decalcify bone without destroying antigen sites, it works especially well to decalcify tissue for immunohistochemistry.

Do not breathe dusts, mists, or vapors. Wash body thoroughly after handling. Wear protective gloves and clothing. Wear eye and face protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off all contaminated clothing. Rinse skin with water/shower. Take off immediately all contaminated clothing and wash it before reuse. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a doctor. Specific treatment, see first aid section on SDS. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing Call a doctor if you feel unwell. STORAGE: Store locked up. DISPOSAL: Dispose of contents/container in accordance with local regulations. NOTICE: Read safety data sheet for additional information.

**ANTIDOTE**

**External:** Flood with water for at least 15 minutes.

**Internal:** Raw egg white, milk, rice gruel. Call a physician.

**DANGER:** Contains dilute acid. May be corrosive to metallic parts, skin, eyes and clothing. Avoid prolonged contact.

EC

REP

mdi Europa GmbH Langenhagener Str. 71 D-30855 Langenhagen Germany +49-511-3908 9530



Cancer Diagnostics, Inc. 4300 Emperor Blvd., 400 Durham, NC 27703 (877) 846 5393 Fax: (877) 817 1716

REF



CF038 Gal (3.8L)



CF100 Liter (1000mL)



CF438 Gal (3.8L), CS/4



CF400 Liter (1000mL), CS/4

