

Lubricant is specially formulated for cleaning, lubricating and protecting microtomes and cryostat machines in high volume laboratories and extreme temperature environments.

#### DIRECTIONS

Insert extension tube. Hold can in upright position. Spray in short even bursts. Can be used with machines turned on or turned off. May be used on all moving parts. Make sure to remove blade before using.

#### CAUTION

FLAMMABLE. CONTENTS UNDER PRESSURE. CONTENTS OF THIS CONTAINER CAN IGNITE UNDER CERTAIN CIRCUMSTANCES. DO NOT USE NEAR FIRE, HEAT, SPARKS, OPEN FLAME OR WHILE SMOKING. DO NOT EXPOSE TO HEAT OR STORE AT TEMPERATURES ABOVE 120 F/49 C. DO NOT STORE CAN IN ENCLOSED VEHICLES, DIRECT SUNLIGHT OR OTHER SOURCES OF HEAT. DO NOT PUNCTURE OR INCINERATE CONTAINER EVEN WHEN EMPTY. USE ONLY WITH ADEQUATE VENTILATION. HARMFUL OR FATAL IF SWALLOWED.

**FIRST AID: Eye Contact:** Flush with water for 15 minutes. If irritation persists call physician. **Skin Contact:** Wash thoroughly with soap and water. If irritation persists call physician. **Inhalation:** Remove to fresh air. If breathing is difficult seek medical attention. **Ingestion:** Do not induce vomiting. May be fatal if swallowed. Call physician immediately.



Flammable Gas



Gas Under Pressure



Acute Toxicity



24 Hour Emergency Contact: **CHEMTREC**, 1-800-424-9300 USA and Canada

Order Toll Free: (877) 846 5393 Fax: (877) 817 1716 [www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)

# Lüben™

## (2 in 1) Cryostat & Microtome Lubricant

For Microtomes & Cryostats  
Low Temperature -55°F

#### KEEP OUT OF REACH OF CHILDREN

CONTENTS UNDER PRESSURE.  
READ BACK OF LABEL CAREFULLY.

#### FOR LABORATORY USE ONLY

3 oz (85g) MCL1000 (ea.)

For in vitro diagnostic use only

Spray Lubricant



CANCER  
DIAGNOSTICS, INC.

Made in USA