



For professional use only.

Flammable

- Intended Use -

Insert extension tube. Place finger on trigger when removing tab from trigger assembly. Tear tab from trigger assembly. 60 below™ rapidly cools paraffin blocks and tissue. Two quick bursts ensure deeper penetration into specimens. 60 below™ contains no CFC's and leaves no residue. Spray 6-8 inches towards surface.

Read entire label before using. Hold can nozzle approximately 3" from paraffin block and tissue. Spray short blasts directly onto surface. Frost will appear on surface when cooled. Hold can upright when spraying. 60 below™ contains a liquefied propellant which will expel as a non-cooling gas if can is held upside down, sideways, or is shaken. If sprayed for extended time, container will become very cold and the pressure will weaken. Spray short blasts or wear gloves to insulate hands.

- Packaging -

Catalog#	Volume
06212	CS/12
06206	CS/6

Non-Flammable

- Intended Use -

Insert extension tube. Place finger on trigger when removing tab from trigger assembly. Tear tab from trigger assembly. 60 below™ rapidly cools paraffin blocks and tissue. Two quick bursts ensure deeper penetration into specimens. 60 below™ contains no CFC's, is nonflammable, and leaves no residue. Spray 6-8 inches towards surface.

Read entire label before using. Hold can nozzle approximately 3" from paraffin block and tissue. Spray short blasts directly onto surface. Frost will appear on surface when cooled. Hold can upright when spraying. 60 below™ contains a liquefied propellant which will expel as a non-cooling gas if can is held upside down, sideways, or is shaken. If sprayed for extended time, container will become very cold and the pressure will weaken. Spray short blasts or wear gloves to insulate hands.

- Packaging -

Catalog#	Volume
06212-NF	CS/12

