

TriCide 10™



Cleaner • Disinfectant • Deodorizer Fungicide • Mildewstat • Virucide*

TriCide 10™ when used as directed is formulated to disinfect, clean, sanitize, & deodorize and is formulated for use on washable hard, non-porous surfaces such as: counters, floors, walls, glass surfaces, painted surfaces, plastic, metal, stainless steel, & medical equipment surfaces. For use in laboratories, hospitals, clean rooms, autopsy rooms, morgues & veterinary clinics.

ACTIVE INGREDIENTS:

| | |
|--|----------|
| Alkyl (50%C14, 40%C12, 10%C16) | |
| dimethyl benzyl ammonium chloride..... | 0.034% |
| Octyl decyl dimethyl ammonium chloride | 0.026% |
| Didecyl dimethyl ammonium chloride..... | 0.013% |
| Dioctyl dimethyl ammonium chloride..... | 0.013% |
| OTHER INGREDIENTS: | 99.914% |
| TOTAL: | 100.000% |

KEEP OUT OF REACH OF CHILDREN CAUTION

See side panel for additional precautionary statements and first aid.

EPA REG. NO. 10324-85-93437 • EPA EST. NO. 5747-MI-1

NET CONTENTS: 1 Quart (946 mL)

DISTRIBUTED BY:



4300 Emperor Blvd., # 400, Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com

DISINFECTANT PERFORMANCE: This product kills the following bacteria in 10 minutes on hard, non-porous surfaces: *Burkholderia cepacia*, *Campylobacter jejuni*, *Corynebacterium ammoniagenes*, *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterobacteriaceae*, *Enterococcus faecalis*, *Enterococcus faecium* (VRE), *Escherichia coli*, *Escherichia coli* O157:H7, *Klebsiella pneumoniae*, *Legionella pneumophila*, *Listeria monocytogenes*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Salmonella enterica*, *Salmonella schottmuelleri*, *Salmonella typhi*, *Serratia marcescens*, *Shigella dysenteriae*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus aureus* (MRSA), *Staphylococcus aureus* (CA MRSA), *Staphylococcus epidermidis*, *Streptococcus pyogenes*, *Vibrio cholerae*.

VIROCIDAL* PERFORMANCE: This product kills the following viruses in 10 minutes on hard, non-porous surfaces: Avian Influenza A (H5N1) Virus, Avian Influenza A/Turkey/Wisconsin Virus, Hepatitis B Virus (HBV), Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus), Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2, Human Coronavirus, Human Immunodeficiency Virus Type 1† (HIV-1), Influenza A (H1N1) Virus, Influenza A/Brazil Virus, Respiratory Syncytial Virus, Vaccinia Virus.

†Indicates a 2-minute contact time is required for this claim.

ANIMAL PREMISE VIROCIDAL* PERFORMANCE: This product kills the following viruses in 10 minutes on hard, non-porous surfaces: Canine Coronavirus, Canine Distemper Virus, Infectious Bovine Rhinotracheitis Virus, Newcastle Disease Virus, Porcine Respiratory & Reproductive Virus, Porcine Rotavirus, Pseudorabies Virus, Transmissible Gastroenteritis Virus.

FUNGICIDAL PERFORMANCE: This product kills the following fungi in 10 minutes on hard, non-porous surfaces: *Candida albicans*, *Trichophyton mentagrophytes*.

MILDEWSTATIC PERFORMANCE: This product controls the following mold on hard, non-porous surfaces: *Aspergillus niger*.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Rinse hard, non-porous food contact surfaces with potable water after application of product. This product is not for use on medical device surfaces.

FOR USE AS A HOSPITAL DISINFECTANT, VIROCIDE*, DEODORIZER, CLEANER:

1. Pre-clean heavily soiled areas.
2. Apply this product undiluted to disinfect hard, non-porous surfaces with a sponge, brush, cloth, hand pump trigger spray device, coarse trigger spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray.
3. Treated surfaces must remain wet for 10 minutes.
4. Wipe dry with a clean cloth or allow to air dry.
5. Change solution daily or when visibly dirty.

***KILLS HIV, HBV AND HCV ON PRE-CLEANED HARD, NON-POROUS SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV AND HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

Personal Protection: Wear protective latex gloves, gowns, masks and eye protection.
Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from hard, non-porous surfaces and objects before application of this product.

Disposal of Infectious Materials: Blood and other body fluids, cleaning materials and clothing must be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

Contact Time: Allow hard, non-porous surface to remain wet for 2 minutes to kill HIV and for 10 minutes to kill all other viruses and bacteria listed on the label.

TO KILL FUNGI: Pre-clean heavily soiled surfaces. Apply product to hard, non-porous surfaces. Allow surface to remain wet for 10 minutes. Wipe surfaces (and let air dry).

TO CONTROL MOLD/MILDEW: Apply product to hard, non-porous surfaces which will effectively inhibit the growth of mold and mildew and their odors. Repeat treatment every seven days, or more often if new growth appears.

ANIMAL PREMISES: Prior to use of this product, remove all animals and feeds from areas to be treated, animal transportation vehicles, and enclosures. Remove all litter, droppings and manure from floors, walls and surfaces of barns, pens, stalls, chutes and other surfaces of facilities and fixtures occupied or traversed by animals. Empty all troughs, racks and other feeding and watering appliances. Thoroughly clean surfaces with soap or detergent and rinse with water.

FOR USE AS AN ANIMAL PREMISE DISINFECTANT/VIROCIDE*: Apply product to disinfect hard, non-porous surfaces with a sponge, brush, cloth, hand pump trigger spray device, coarse trigger spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray. Immerse all halters and other types of equipment used in handling and restraining animals, as well as forks, shovels, and scrapers used for removing litter and manure in the use solution. Treated surfaces must remain wet for 10 minutes. Ventilate buildings, coops and other closed spaces. Do not house animals or employ equipment until treatment has been absorbed, set or dried. Thoroughly scrub all treated feed racks, troughs, automatic feeders, fountains and waterers and other treated equipment which can contact food or water with soap or detergent, and rinse with potable water before reuse.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

STORAGE: Store only in original container. Keep this product under locked storage sufficient to make it inaccessible to children or persons unfamiliar with its proper use.

PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.

CONTAINER HANDLING: Non-Refillable Container. Do not reuse or refill this container. Triple rinse container promptly after emptying. Triple rinse as follows: Fill the container ¼ full with water and recap. Shake for 10 seconds. Drain for 10 seconds after the flow begins to drip. Follow Pesticide Disposal instructions for rinsate disposal. Repeat this procedure two more times. Then offer for recycling if available or puncture and dispose of in a sanitary landfill, or by incineration.

PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION. Causes moderate eye irritation. Harmful if swallowed or absorbed through the skin. Avoid contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

ENVIRONMENTAL HAZARDS

This product is toxic to fish and freshwater aquatic invertebrates.

FIRST AID

In case of emergency, call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

IF ON SKIN: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.

IF SWALLOWED: Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give anything by mouth to an unconscious person. Call a poison control center or doctor immediately for treatment advice.



TriCide 10TM

DISTRIBUTED BY:



CANCER
DIAGNOSTICS, INC.

4300 Emperor Blvd., # 400, Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

EPA REG. NO. 10324-85-93437 • EPA EST. NO. 5747-MI-1

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION. Causes moderate eye irritation. Harmful if swallowed or absorbed through the skin. Avoid contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

ENVIRONMENTAL HAZARDS

This product is toxic to fish and freshwater aquatic invertebrates.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

PESTICIDE STORAGE: Store only in original container. Keep this product under locked storage sufficient to make it inaccessible to children or persons unfamiliar with its proper use.

PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.

CONTAINER HANDLING: Non-Refillable Container. Do not reuse or refill this container. Triple rinse container (or equivalent) promptly after emptying. Triple rinse as follows: Fill the container ¼ full with water and recap. Shake for 10 seconds. Drain for 10 seconds after the flow begins to drip. Follow Pesticide Disposal instructions for rinsate disposal. Repeat this procedure two more times. Then offer for recycling if available or puncture and dispose of in a sanitary landfill, or by incineration.

ACTIVE INGREDIENTS:

| | |
|--|-----------------|
| Alkyl (50%C14, 40%C12, 10%C16) dimethyl benzyl ammonium chloride..... | 0.034% |
| Octyl decyl dimethyl ammonium chloride | 0.026% |
| Didecyl dimethyl ammonium chloride | 0.013% |
| Diocetyl dimethyl ammonium chloride..... | 0.013% |
| OTHER INGREDIENTS: | 99.914% |
| TOTAL: | 100.000% |

KEEP OUT OF REACH OF CHILDREN CAUTION

See side panel for additional precautionary statements.

FIRST AID

In case of emergency, call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

IF ON SKIN: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.

IF SWALLOWED: Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give anything by mouth to an unconscious person. Call a poison control center or doctor immediately for treatment advice.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

See inside for other directions for use and organisms.



NET CONTENTS: 12X1 QUART (TOTAL: 3 GALLONS)