

**\*VIRUCIDAL ACTIVITY:** This product is an effective virucide on hard, non-porous surfaces against: Human Immunodeficiency Virus Type 1 (HIV-1) when the treated surface is allowed to remain wet for 2 minutes. Norovirus (Norwalk Virus) and Feline calicivirus when the treated surface is allowed to remain wet for 10 minutes. Influenza A (H1N1) Virus when the treated surface is allowed to remain wet for 10 minutes. Use additional wipes if needed to assure that the treated surface remains continuously wet for the prescribed time.

**KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.**

**PERSONAL PROTECTION:** When handling items soiled with blood or body fluids use disposable latex gloves, gowns, masks, and eye coverings.

**CLEANING PROCEDURES:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this product. Wipes can be used for this purpose.

**CONTACT TIME:** Thoroughly wet surface with a new wipe. Allow surface to air dry. Although efficacy at a 2 minute contact time has been shown to be adequate against HIV-1, use the contact time for the disinfection of other organisms listed on this label.

**DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

Not for personal cleansing. • This is not a baby wipe! • Do not flush in toilet.

**TO OPEN:**

Remove lid from canister. Lift flap in center of lid. Pull up corner of center wipe, twist to a point and thread through dispenser hole in lid. DO NOT PUSH FINGER THROUGH OPENING. Replace lid. Pull out exposed wipe and snap off. The next wipe pops up automatically. Dispose of used wipes in trash receptacle. Snap center flap down when finished to retain moisture. Do not flush.

**PRECAUTIONARY STATEMENTS**

Hazards to Humans and Domestic Animals

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum or using tobacco.

**FIRST AID**

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

**IF IN EYES:** Hold eyes open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.


**STORAGE AND DISPOSAL**

Do not contaminate water, food or feed by storage or disposal.

**PESTICIDE STORAGE:** Store in an area inaccessible to children.

**CONTAINER STORAGE/DISPOSAL:** Wipes: Dispose of wipe in trash after use. Do not flush. Canister: Tightly close lid between uses to retain moisture. Nonrefillable container. Do not reuse or refill this container. Offer empty container for recycling if available or discard in trash.

# TriCide<sup>10</sup>™ Wipes



**\*VIRUCIDAL · TUBERCULOCIDAL  
BACTERICIDAL · PSEUDOMONICIDAL**

**ACTIVE INGREDIENTS**

Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> ) dimethyl benzyl ammonium chlorides . . . . .	0.105%
Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides . . . . .	0.105%
<b>INERT INGREDIENTS</b> . . . . .	<u>99.790%</u>
<b>Total</b> . . . . .	<b>100.000%</b>

## KEEP OUT OF REACH OF CHILDREN CAUTION

See side panels for additional precautionary statements.

EPA REG. NO. 1839-174-93437 · EPA EST. NO. 5747-MI-01

NET CONTENTS: 2 lbs. 5 oz. · 160 Pre-measured wipes 6" X 6 ¾"

U.S. Patent No. 5,444,094

This product is bactericidal in the presence of organic soil (5% blood serum) against Escherichia coli (E. coli), Escherichia coli O157:H7, Salmonella (choleraesuis) enterica, Pseudomonas aeruginosa, Vancomycin Resistant Enterococcus faecalis (VRE), Staphylococcus aureus and Methicillin Resistant Staphylococcus aureus (MRSA) at 10 minutes contact time.

This product is bactericidal in the presence of organic soil (5% blood serum) against Vancomycin Intermediate Resistant Staphylococcus aureus (VISA), Clinical Isolate No. 99286, at 10 minutes contact time.

This product is tuberculocidal in the presence of organic soil (5% blood serum) against Mycobacterium bovis BCG (Tuberculosis) at 20°C at 10 minutes contact time.

This product is virucidal in the presence of organic soil (5% blood serum) against Human Immunodeficiency Virus Type 1, Strain HTLV-IIIb (HIV-1) at 2 minutes contact time and against Norovirus (Norwalk Virus) and Feline Calicivirus (FCV) at 10 minutes contact time.

This product is virucidal in the presence of organic soil (5% blood serum) against Norovirus (Norwalk Virus) and Feline Calicivirus (FCV) at 10 minutes contact time.

This product is a non-woven, disposable wipe, pre-saturated with a quaternary disinfectant cleaner.

This product is intended for use in: hospitals and healthcare settings: ambulances, anesthesia, CAT labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's offices, Emergency Rooms, Home health care, Hospices, ICU, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating rooms, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physician's offices, Ophthalmic offices, Radiology, Recovery Room, Rehabilitation, Respiratory therapy, Surgery, Transport vehicles, X-Ray, Veterinary clinics, and animal hospitals.

This product will clean and disinfect hard, non-porous, surfaces of: bed railings, cabinets, carts, chairs, counters, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats, toilets, bathtubs, bathroom fixtures, shower stalls, sinks, tiles, vanity tops, hospital beds, traction devices, MRI, CAT, scales, paddles and wheelchairs.

**DIRECTIONS FOR USE**

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is not for use on critical and semi-critical medical device surfaces.

**CLEANING:**

Use one wipe to completely remove surface of all gross debris from surfaces to be treated.

**DISINFECTING:**

To disinfect hard, non-porous surfaces, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for disinfection: To disinfect against Escherichia coli (E. coli), Escherichia coli O157:H7, Pseudomonas aeruginosa, Salmonella (choleraesuis) enterica and Vancomycin Resistant Enterococcus faecalis (VRE), allow the treated surface to remain wet for 10 minutes. To disinfect against Mycobacterium bovis BCG (Tuberculosis), allow the treated surface to remain wet for 10 minutes. To disinfect against Staphylococcus aureus and Methicillin Resistant Staphylococcus aureus (MRSA), allow the treated surface to remain wet for 10 minutes. To disinfect against Vancomycin intermediate Resistant Staphylococcus aureus (VISA), allow the treated surface to remain wet for 10 minutes. A potable water rinse is required for surfaces in direct contact with food. For plastic and painted surfaces, spot test on an inconspicuous area before use. For heavily soiled surfaces, clean before following disinfecting instructions.

Distributed by:



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com