



TriCide10™ Wipes

- Intended Use -

KEEP OUT OF REACH OF CHILDREN

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Rinse hard, non-porous food contact surfaces with potable water after application of product. This product is not for use on critical and semi-critical medical device surfaces. Not for personal cleansing. This is **NOT** a baby wipe! Do not flush in toilet.

- Active Ingredients -

Alkyl (60% C ₁₄ , 30% C ₁₆ , 5% C ₁₂ , 5% C ₁₈) dimethyl benzyl ammonium chlorides	0.105%
Alkyl (68% C ₁₂ , 32% C ₁₄) dimethyl ethylbenzyl ammonium chlorides	0.105%
Inert Ingredients	99.790%
TOTAL	100.000%

- Cleaning & Disinfecting -

To clean, use one wipe to completely remove surface of all gross debris from surfaces to be treated. To disinfect hard, non-porous surfaces, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for disinfection: To disinfect against Escherichia coli (E. coli), Escherichia coli O157:H7, Pseudomonas aeruginosa, Salmonella (choleraesuis) enterica and Vancomycin Resistant Enterococcus faecalis (VRE), allow the treated surface to remain wet for 10 minutes. To disinfect against Mycobacterium bovis BCG (Tuberculosis), allow the treated surface to remain wet for 10 minutes. To disinfect against Staphylococcus aureus and Methicillin Resistant Staphylococcus aureus (MRSA), allow the treated surface to remain wet for 10 minutes. To disinfect against Vancomycin intermediate Resistant Staphylococcus aureus (VISA), allow the treated surface to remain wet for 10 minutes. A potable water rinse is required for surfaces in direct contact with food. For plastic and painted surfaces, spot test on an inconspicuous area before use. For heavily soiled surfaces, clean before following disinfecting instructions.

- Special Instructions -

For cleaning and decontamination against HIV-1 on surfaces/objects soiled with blood/body fluids.

Personal Protection	When handling items soiled with blood or body fluids use disposable latex gloves, gowns, masks, and eye coverings.
Cleaning Procedure	Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this product. Wipes can be used for this purpose.
Disposal of Infectious Materials	Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.
Contact Time	Thoroughly wet surface with a new wipe. Allow surface to air dry. Although efficacy at a 2 minute contact time has been shown to be adequate against HIV-1, use the contact time for the disinfection of other organisms listed on product label.

- To Open -

Remove lid from canister. Lift flap in center of lid. Pull up corner of center wipe, twist to a point and thread through dispenser hole in lid. DO NOT PUSH FINGER THROUGH OPENING. Replace lid. Pull out exposed wipe and snap off. The next wipe pops up automatically. Dispose of used wipes in trash receptacle. Snap center flap down when finished to retain moisture. Do not flush.

- Packaging -

Catalog#	Volume	Qty.
TC10W-4	Canister ^(160 wipes)	CS/4

