




Luxol Fast Blue Stain Kit

Description: The Luxol Fast Blue Stain Kit is designed for staining myelin/myelinated axons and Nissl substance on formalin fixed, paraffin-embedded tissue as well as frozen tissue. This product is used for identifying the basic neuronal structure in brain or spinal cord sections.

Myelinated Fibers: Blue
Nissl Substance: Violet
Nerve Cells: Violet

Uses/Limitations: Not to be taken internally.
For In-Vitro Diagnostic use only.
Histological applications.
Do not use if reagents become cloudy.
Do not use past expiration date.
Use caution when handling reagents.
Non-Sterile.


Control Tissue: Cerebral Cortex
Spinal Cord



Kit Contents:

| <u>Kit Contents</u> | <u>Volume</u> | <u>Storage</u> |
|------------------------------------|---------------|----------------|
| Cresyl Echt Violet Solution | 125 ml | 2-8°C |
| Luxol Fast Blue Solution | 125 ml | 18-25°C |
| Lithium Carbonate Solution (0.05%) | 500 ml | 18-25°C |
| Alcohol, Reagent (70%) | 500 ml | 18-25°C |

Precautions: Avoid contact with skin and eyes.
Harmful if swallowed.
Follow all Federal, State, and local regulations regarding disposal.

Storage: 2° C  25° C

**Mixed Storage Conditions.
Separate Contents.**




**Procedure:**

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Incubate slide in Luxol Fast Blue Solution for 24 hours at room temperature or 2 hours at 60°C.
3. Rinse thoroughly in distilled water.
4. Differentiate section by dipping in Lithium Carbonate Solution (0.05%) several times (up to 20 seconds).
5. Continue differentiation by repeatedly dipping in Alcohol, Reagent (70%) until gray-matter is colorless and white-matter remains blue.
6. Rinse slide in 2 changes of distilled water.
7. Incubate slide in Cresyl Echt Violet (0.1%) for 2-5 minutes.
8. Rinse quickly in 1 change of distilled water.
9. Dehydrate quickly in 3 changes of absolute alcohol.
10. Clear as desired and mount in synthetic resin.

References:

1. Sheenan, D.C., Hrapchak, B.B. Theory and Practice of Histotechnology, 2nd Edition. Battelle Press, Columbus, OH. Page 262-264. 1980
2. Kluver, H., Barrera, E.A. A Method for the combined staining of cells and fibers in the nervous system. Journal of Neuropathology and Experimental Neurology, 1953, 12: pages 400-403.

| Description: | Volume |
|------------------------------------|---------|
| Cresyl Echt Violet Solution | 125 ml |
| | 500 ml |
| | 1000 ml |
| Luxol Fast Blue Solution | 125 ml |
| | 500 ml |
| | 1000 ml |
| Lithium Carbonate Solution (0.05%) | 500 ml |
| | 1000 ml |
| Alcohol, Reagent (70%) | 500 ml |
| | 1000 ml |

Storage: 2° C  25° C**Mixed Storage Conditions.
Separate Contents.**