



# SuperLite™ Storage Units

## Tissue Cassette and Microscope Slides



### TISSUE CASSETTES

#### TWO VERSIONS

- 1 Plastic Drawer Inserts
- 2 Corrugated Partitions



### MICROSCOPE SLIDES

#### TWO VERSIONS

- 1 Plastic Drawer Inserts
- 2 Corrugated Partitions



### CASSETTE STORAGE ITEM# CPC-6



### CASSETTE STORAGE ITEM# ACC-6



- 1 12 each, high impact polystyrene plastic trays, (7 1/2" wide x 13 1/4" deep x 1 3/4" high) /or Corrugated Partitioning
- 2 12 each, sturdy corrugated drawers with front labeling for easy retrieval
- 3 6 each, corrugated drawer enclosures for dust-proof protection
- 4 Reinforced Kraft Outer Carton for long term storage
- 5 Case dimensions: 27" wide x 16 1/2" deep x 10" high

### SLIDE STORAGE ITEM# CPS-6



### SLIDE STORAGE ITEM# ACS-6



- 1 6 each, high impact polystyrene plastic trays, (7 1/2" wide x 13 1/4" deep x 2 5/8" high) /or Corrugated Partitioning
- 2 6 each, sturdy corrugated drawers with front labeling for easy retrieval
- 3 6 each, corrugated drawer enclosures for dust-proof protection
- 4 Reinforced Kraft Outer Carton for long term storage
- 5 Case dimensions: 27" wide x 16 1/2" deep x 10" high



**Cancer Diagnostics Inc.**

**Contact Information**

**Visit us online**

4300 Emperor Blvd-400  
Durham, North Carolina 27703

**Toll Free:** (877) 846 5393  
**Fax:** (877) 817 1716

**cancerdiagnostics.com**

©2016 Cancer Diagnostics, Inc. All rights reserved. Printed in the USA