



Glacial Acetic Acid

#C8551 500ML

#C8552 1 GAL (3.8L)

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____



DANGER HAZARD STATEMENTS

Toxic if inhaled. Harmful in contact with skin. Causes severe skin burns and serious eye damage. Flammable liquid and vapor.

PRECAUTIONARY STATEMENTS - PREVENTION

Use only outdoors or in a well-ventilated area. Wear protective gloves and clothing. Do not breathe dusts, vapors, or mists. Wash body thoroughly after handling. Wear eye and face protection. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing.

IF ON SKIN: Wash with plenty of water. Call a doctor if you feel unwell. Take off immediately all contaminated clothing and wash it before reuse. Immediately call a doctor.

IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a doctor. Specific treatment, see first aid section on SDS.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store in a well-ventilated place. Keep container tightly closed. Store locked up. Keep cool. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA