

Revision Date 2024-07-26

Version 1

## 1. IDENTIFICATION

**Product Name** Michel Transport Medium  
**Product Code** FX1046, FX1046M, FX1046-500, FX1046-20, FX1046-96  
**Recommended Use** For laboratory, scientific, R&D or manufacturing use.  
**Company** Cancer Diagnostics, Inc.  
116 Page Point Circle  
Durham, NC 27703  
Tel. (877) 846-5393  
**Emergency Telephone** Call CHEMTREC 1-800-424-9300

## 2. HAZARDS IDENTIFICATION

### Classification

#### **OSHA Regulatory Status**

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Acute toxicity - Oral

Category 4

### Label elements

#### **Signal word**

Warning

#### **Hazard statements**

Harmful if swallowed.



#### **Precautionary Statements - Prevention**

Wash face, hands and any exposed skin thoroughly after handling. Do not eat, drink or smoke when using this product.

IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Rinse mouth.

#### **Precautionary Statements - Storage**

Store tightly closed at 15-30C.

#### **Precautionary Statements - Disposal**

Dispose of contents/container to an approved waste disposal plant.

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS No.	Weight-%
Water	7732-18-5	80-90
Ammonium sulfate	7783-20-2	35-45
Potassium citrate	6100-05-6	<2
n-Ethylmaleimide	128-53-0	<2
Magnesium sulfate, heptahydrate	10034-99-8	<2

### 4. FIRST AID MEASURES

#### Description of first aid measures

<b>Eye contact</b>	Flush eyes with plenty of water, removing contact lenses if present. Get medical attention if irritation develops.
<b>Skin contact</b>	Wash thoroughly with soap and water while removing contaminated garments. Get medical attention if irritation develops.
<b>Inhalation</b>	Remove to fresh air. Get medical attention for any breathing difficulty.
<b>Ingestion</b>	Rinse mouth and drink several glasses of water. Contact a physician or poison control center if symptoms develop.

#### Most important symptoms and effects, both acute and delayed

**Symptoms** No information available.

### 5. FIRE-FIGHTING MEASURES

#### Suitable extinguishing media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment

#### Specific hazards arising from the chemical

No information available.

#### Hazardous combustion products

Carbon oxides. Ammonia. Nitrogen oxides (NOx). Sulphur oxides.

#### Protective equipment and precautions for firefighters

Firefighters should wear self-contained breathing apparatus with full facepiece operated in pressure-demand or other positive pressure mode.

**NFPA**                      **Health hazards** 1                      **Flammability** 0                      **Instability** 0                      **Physical and Chemical Properties** -

### 6. ACCIDENTAL RELEASE MEASURES

#### Personal precautions, protective equipment and emergency procedures

**Personal precautions** Use personal protective equipment as required. Ensure adequate ventilation, especially in confined areas.

**Environmental precautions** Prevent entry into waterways, sewers, basements or confined areas.

#### Methods and material for containment and cleaning up

**Methods for containment** Prevent further leakage or spillage if safe to do so.

**Methods for cleaning up** Pick up and transfer to properly labeled containers.

## 7. HANDLING AND STORAGE

**Precautions for safe handling** Use personal protective equipment as required  
Handle in accordance with good industrial hygiene and safety practice.

**Storage Conditions** Keep containers tightly closed in a dry, cool and well-ventilated place.

**Incompatible materials** Strong acids. Strong bases. Oxidizers.

## 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

**Occupational exposure limits** This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies

### Appropriate engineering controls

**Engineering Controls** Emergency showers, eyewash stations, ventilation systems.

### Individual protection measures, such as personal protective equipment

**Eye/face protection** Wear safety glasses with side shields (or goggles).

**Skin and body protection** Wear protective gloves and protective clothing.

**Respiratory protection** If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

**General Hygiene Considerations** Handle in accordance with good industrial hygiene and safety practice.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

### Information on basic physical and chemical properties

<b>Physical state</b>	Liquid
<b>Appearance</b>	Clear, colorless - yellow
<b>Odor</b>	Pungent
<b>Odor threshold</b>	No information available
<b>pH</b>	7.0 - 7.2
<b>Melting point / freezing point</b>	No information available
<b>Boiling point / boiling range</b>	No information available
<b>Flash point</b>	No information available
<b>Evaporation rate</b>	No information available
<b>Flammability (solid, gas)</b>	No information available
<b>Flammability Limit in Air</b>	
<b>Upper flammability limit:</b>	No information available
<b>Lower flammability limit:</b>	No information available
<b>Vapor pressure</b>	No information available
<b>Vapor density</b>	No information available
<b>Relative density</b>	No information available
<b>Water solubility</b>	Miscible with water
<b>Solubility in other solvents</b>	No information available
<b>Partition coefficient</b>	No information available
<b>Autoignition temperature</b>	No information available
<b>Decomposition temperature</b>	No information available
<b>Kinematic viscosity</b>	No information available

## 10. STABILITY AND REACTIVITY

<b>Stability</b>	Stable under recommended storage conditions.
<b>Possibility of Hazardous Reactions</b>	None under normal processing.
<b>Conditions to avoid</b>	Extremes of temperature and direct sunlight.
<b>Incompatible materials</b>	Strong acids. Strong bases. Oxidizers.
<b>Hazardous Decomposition Products</b>	Carbon oxides. Nitrogen oxides (NOx). Sulphur oxides. Ammonia.

## 11. TOXICOLOGICAL INFORMATION

### Information on likely routes of exposure

<b>Inhalation</b>	No data available.
<b>Eye contact</b>	No data available.
<b>Skin contact</b>	No data available.
<b>Ingestion</b>	No data available.

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Ammonium sulfate 7783-20-2	= 2000 mg/kg ( Rat )	-	-
n-Ethylmaleimide 128-53-0	= 25 mg/kg ( Rat )	-	-

### Delayed and immediate effects as well as chronic effects from short and long-term exposure

<b>Sensitization</b>	No information available.
<b>Germ cell mutagenicity</b>	No information available.
<b>Carcinogenicity</b>	This chemical does not contain any carcinogens or potential carcinogens as listed by ACGIH, OSHA, IARC or NTP

## 12. ECOLOGICAL INFORMATION

### Ecotoxicity

Chemical Name	Algae/aquatic plants	Fish	Crustacea
Ammonium sulfate 7783-20-2	-	250: 96 h Brachydanio rerio mg/L LC50 480: 96 h Brachydanio rerio mg/L LC50 flow-through 420: 96 h Brachydanio rerio mg/L LC50 semi-static 18: 96 h Cyprinus carpio mg/L LC50 32.2 - 41.9: 96 h Oncorhynchus mykiss mg/L LC50 flow-through 5.2 - 8.2: 96 h Oncorhynchus mykiss mg/L LC50 static 100: 96 h Pimephales promelas mg/L LC50 123 - 128: 96 h Poecilia reticulata mg/L LC50 semi-static 126: 96 h Poecilia reticulata mg/L LC50 460 - 1000: 96 h Leuciscus idus mg/L LC50 static	14: 48 h Daphnia magna mg/L LC50 423: 24 h Daphnia magna mg/L EC50

### Persistence and degradability

No information available.

**Bioaccumulation**

No information available.

Chemical Name	Partition coefficient
Ammonium sulfate 7783-20-2	-5.1

**Other adverse effects** No information available

**13. DISPOSAL CONSIDERATIONS**

**Waste treatment methods**

**Disposal of wastes** Disposal should be in accordance with applicable regional, national and local laws and regulations.

**Contaminated packaging** Do not reuse container. Emptied containers may contain residue. Continue to follow label warnings after container is emptied.

**14. TRANSPORT INFORMATION**

Transportation information is provided as a general reference only and may not be applicable in all situations. This information applies to non-bulk shipments only. Per 49 CFR §173.22, it is the shipper's responsibility to ensure that all materials are properly packaged, classified and labeled prior to shipment.

**DOT** Not regulated

**IATA** Not regulated

**15. REGULATORY INFORMATION**

**US Federal Regulations**

**SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product contains a chemical or chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

Chemical Name	SARA 313 - Threshold Values %
Ammonium sulfate - 7783-20-2	1.0

**SARA 311/312 Hazard Categories**

<b>Acute health hazard</b>	Yes
<b>Chronic Health Hazard</b>	No
<b>Fire hazard</b>	No
<b>Sudden release of pressure hazard</b>	No
<b>Reactive Hazard</b>	No

**CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

**CERCLA**

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level

pertaining to releases of this material

**US State Regulations**

**California Proposition 65**

This product does not contain any Proposition 65 chemicals

**U.S. State Right-to-Know Regulations**

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Ammonium sulfate 7783-20-2	-	X	X

**16. OTHER INFORMATION**

**Prepared By**  
**Revision Date**

CDI Regulatory Affairs (Email: [compliance@eki-chem.com](mailto:compliance@eki-chem.com))  
2024-07-26

**Disclaimer**

7 UbWf 8 ]U[ bcghWz-bW makes no representation as to the comprehensiveness or accuracy of this document. Individuals using this information, or the product to which it refers, must exercise their independent judgment in determining all appropriateness for a particular purpose. Accordingly, 7 UbWf 8 ]U[ bcghWz-bW will not be responsible for damages of any kind or nature resulting from the use of this information or the corresponding product. No representations or warranties of any kind or nature, including but not limited to: express warranties, implied warranties or merchantability, or warranties of fitness for a particular purpose, are made hereunder with respect to the information set forth herein or to the product to which the information refers.

**End of Safety Data Sheet**