



Bile Stain Kit, Hall

Solutions Provided:

#SSC1076-250 - **Fouchet's Reagent**

#SSC1135-250 - **Van Gieson's Stain**

Conventional Procedure

1. Deparaffinize and hydrate sections through alcohol to distilled water.
2. Rinse slides thoroughly in distilled water.
3. Stain sections in Fouchet's Reagent for 5 minutes.
4. Rinse in 2-3 changes of tap or distilled water.
5. Stain sections in Van Gieson's Stain for 5 minutes.
6. Wash in running tap water for 5 minutes.
7. Rinse slides in 1 change of 95% alcohol.
8. Dehydrate slides in 2-3 changes of 100% alcohol.
9. Clear in 2-3 changes of xylene and mount with appropriate mounting medium.

Results:

Bile or Bilirubin:

Green

Background:

Yellow

