

SAFETY DATA SHEET

Revision Date 2021-08-13

Version 7

1. IDENTIFICATION			
Product Name	10% Formalin Fixative, Neutral, Phosphate Buffered		
Product Code	Formalin, NBF, 10% NBF, (See Addendum for full code listing)		
Recommended Use	Routine histology fixative. For laboratory and in-vitro diagnostic use only.		
Company	Cancer Diagnostics 116 Page Point Circle Durham, NC 27703 Tel. (877) 846-5393		
Emergency Telephone	Call CHEMTREC 1-800-424-9300		

2. HAZARDS IDENTIFICATION

Classification

OSHA Regulatory Status

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Acute toxicity - Inhalation (Gases)	Category 4
Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 1
Skin sensitization	Category 1
Carcinogenicity	Category 1A
Specific target organ toxicity (single exposure)	Category 1
Flammable liquids	Category 4

Label elements

Signal word

Danger

Hazard statements

Harmful if inhaled. Causes skin irritation. Causes serious eye damage. May cause an allergic skin reaction. May cause cancer. Causes damage to organs. Combustible liquid.



Precautionary Statements - Prevention

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Use personal protective equipment as required. Use only outdoors or in a well-ventilated area. Wash face, hands and any exposed skin thoroughly after handling. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves. Do not

breathe dust/fume/gas/mist/vapors/spray. Do not eat, drink or smoke when using this product. Keep away from heat/sparks/open flames/hot surfaces. — No smoking.

Precautionary Statements - Response

IF exposed: Call a POISON CENTER or doctor/physician.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

IF ON SKIN: Wash with plenty of soap and water. Take off contaminated clothing and wash before reuse. If skin irritation or rash occurs: Get medical advice/attention.

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or doctor/physician if you feel unwell.

In case of fire: Use CO2, dry chemical, or foam for extinction.

Precautionary Statements - Storage

Store locked up. Store in a well-ventilated place. Keep cool.

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS No.	Weight-%
Water	7732-18-5	>90
Formaldehyde	50-00-0	3.5-4.0
Methyl alcohol	67-56-1	1.0-1.5
Sodium Phosphate, dibasic	7558-79-4	<1
Sodium Phosphate, monobasic	7558-80-7	<1

4. FIRST AID MEASURES

Description of first aid measures

Eye contact	Immediately flush with plenty of water for at least 15 minutes, separating eyelids occasionally. Remove contact lenses if present. Get immediate medical attention.		
Skin contact	Wash thoroughly with soap and water while removing contaminated garments. Get medical attention if irritation develops.		
Inhalation	Remove to fresh air. Get medical attention for any breathing difficulty.		
Ingestion	Do NOT induce vomiting unless instructed to do so by medical personnel. If conscious, rinse mouth and give several glasses of water to drink. Never give anything by mouth to an unconscious person. Get immediate medical attention.		
Most important symptoms and effects, both acute and delayed			

Symptoms Causes skin, eye and respiratory tract irritation. If swallowed, causes GI disturbances. May cause allergic skin reaction.

5. FIRE-FIGHTING MEASURES

Suitable extinguishing media

CO2, dry chemical, dry sand, alcohol-resistant foam Use extinguishing measures that are appropriate to local circumstances and the surrounding environment

Specific hazards arising from the chemical

May cause sensitization by skin contact. The product causes irritation of eyes, skin and mucous membranes.

Hazardous combustion products

Carbon oxides. Formaldehyde.

Protective equipment and precautions for firefighters

Firefighters should wear self-contained breathing apparatus with full facepiece operated in pressure-demand or other positive pressure mode.

<u>NFPA_</u>	Health hazards 2	Flammability 1	Instability 0	Physical and Chemical Properties -		
6. ACCIDENTAL RELEASE MEASURES						
Personal precautions, protective equipment and emergency procedures						
Personal precaution			othing. Ensure adequate ective equipment as req	e ventilation, especially in uired.		
Environmental prec		Prevent product from entering drains. Should not be released into the environment. Prevent entry into waterways, sewers, basements or confined areas.				
Methods and mater	ial for containment and clea	aning up				
Methods for contain	nment Prevent fu	Prevent further leakage or spillage if safe to do so.				
Methods for cleaning	ng up Absorb spi	Absorb spill with inert material, scoop up and containerize for disposal.				
	7. HANDLING AND STORAGE					
Precautions for safe	•	nal protective equipmen accordance with good ir	t as required idustrial hygiene and sat	fety practice.		
Storage Conditions	Keep cont	ainer tightly closed in a	dry and well-ventilated p	lace. Store at 15C to 25C.		
Incompatible mater	ials Oxidizers.	Oxidizers. Bases.				

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Occupational exposure limits

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Formaldehyde	Ceiling: 0.3 ppm	TWA: 0.75 ppm	IDLH: 20 ppm
50-00-0	0 11	(vacated) TWA: 3 ppm unless	Ceiling: 0.1 ppm 15 min
		specified in 1910.1048	TWA: 0.016 ppm
		(vacated) STEL: 10 ppm 30 min	
		unless specified in 1910.1048	
		(vacated) Ceiling: 5 ppm_unless	
		specified in 1910.1048	
		STEL: 2 ppm see 29 CFR	
		1910.1048	
Methyl alcohol	STEL: 250 ppm	TWA: 200 ppm	IDLH: 6000 ppm
67-56-1	TWA: 200 ppm	TWA: 260 mg/m ³	TWA: 200 ppm
	Skin	(vacated) TWA: 200 ppm	TWA: 260 mg/m ³
		(vacated) TWA: 260 mg/m ³	STEL: 250 ppm
		(vacated) STEL: 250 ppm	STEL: 325 mg/m ³
		(vacated) STEL: 325 mg/m ³	
		(vacated) Skin	

Appropriate engineering controls

Engineering Controls

Emergency showers, eyewash stations, ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/face protection	Wear safety glasses with side shields (or goggles).
Skin and body protection	Wear protective gloves and protective clothing.
Respiratory protection	If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.
General Hygiene Considerations	Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Liquid Clear, colorless Odor of formaldehyde No information available
 7.0 No information available No information available 85 C No information available
No information available Miscible with water
No information available No information available No information available No information available No information available

10. STABILITY AND REACTIVITY

Stability	Stable under recommended storage conditions.
Possibility of Hazardous Reactions	None under normal processing.
Conditions to avoid	Extremes of temperature and direct sunlight.
Incompatible materials	Oxidizers. Bases.
Herendeve Decomposition Dredvets	Formoldobudo, Corbon monovido, Corbon diovido (CO2

Hazardous Decomposition Products Formaldehyde. Carbon monoxide. Carbon dioxide (CO2).

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Inhalation

Harmful by inhalation.

Eye contact	Risk of serious damage to eyes.
-------------	---------------------------------

Skin contact May cause sensitization by skin contact.

Ingestion No data available.

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Formaldehyde 50-00-0	= 500 mg/kg (Rat)	= 270 mg/kg (Rabbit)	= 0.578 mg/L (Rat)4 h
Methyl alcohol 67-56-1	= 5628 mg/kg (Rat)	-	= 83.2 mg/L (Rat)4 h
Sodium Phosphate, dibasic 7558-79-4	= 17 g/kg (Rat)	-	-
Sodium Phosphate, monobasic 7558-80-7	= 8290 mg/kg (Rat)	> 7940 mg/kg (Rabbit)	-

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation Serious eye damage/eye Sensitization Germ cell mutagenicity Carcinogenicity	irritation Risk of serio May cause s No informatio	Irritating to skin. Risk of serious damage to eyes. Irritating to eyes. May cause sensitization by skin contact. No information available. The table below indicates whether each agency has listed any ingredient as a carcinogen		
Chemical Name	ACGIH	IARC	NTP	OSHA
Formaldehyde 50-00-0	A2	Group 1	Known	X
Legend ACGIH (American Conference of Governmental Industrial Hygienists)				

ACGIN (American Conference of Governmental Indus A2 - Suspected Human Carcinogen IARC (International Agency for Research on Cancer) Group 1 - Carcinogenic to Humans NTP (National Toxicology Program)

Known - Known Carcinogen

OSHA (Occupational Safety and Health Administration of the US Department of Labor)

X - Present

12. ECOLOGICAL INFORMATION

Ecotoxicity

Chemical Name	Algae/aquatic plants	Fish	Crustacea
Formaldehyde	-	22.6 - 25.7: 96 h Pimephales	2: 48 h Daphnia magna mg/L LC50
50-00-0		promelas mg/L LC50 flow-through	11.3 - 18: 48 h Daphnia magna
		1510: 96 h Lepomis macrochirus	mg/L EC50 Static
		µg/L LC50 static 41: 96 h	
		Brachydanio rerio mg/L LC50 static	
		0.032 - 0.226: 96 h Oncorhynchus	
		mykiss mL/L LC50 flow-through 100	
		- 136: 96 h Oncorhynchus mykiss	
		mg/L LC50 static 23.2 - 29.7: 96 h	
		Pimephales promelas mg/L LC50	
		static	
Methyl alcohol	-	28200: 96 h Pimephales promelas	-
67-56-1		mg/L LC50 flow-through 100: 96 h	
		Pimephales promelas mg/L LC50	
		static 19500 - 20700: 96 h	
		Oncorhynchus mykiss mg/L LC50	
		flow-through 18 - 20: 96 h	
		Oncorhynchus mykiss mL/L LC50	
		static 13500 - 17600: 96 h Lepomis	
		macrochirus mg/L LC50	
		flow-through	

Persistence and degradability

No information available.

Bioaccumulation

No information available.

Chemical Name	Partition coefficient
Formaldehyde 50-00-0	0.35
Methyl alcohol 67-56-1	-0.77

Other adverse effects

No information available

13. DISPOSAL CONSIDERATIONS

Waste treatment methods

Disposal of wastes Disposal should be in accordance with applicable regional, national and local laws and regulations.

Contaminated packaging

Do not reuse container. Emptied containers may contain residue. Continue to follow label warnings after container is emptied.

Chemical Name	RCRA	RCRA - Basis for Listing	RCRA - D Series Wastes	RCRA - U Series Wastes
Formaldehyde 50-00-0	U122	Included in waste streams: K009, K010, K038, K040, K156, K157	-	U122
Methyl alcohol 67-56-1	-	Included in waste stream: F039	-	U154

Chemical Name	California Hazardous Waste Status
Formaldehyde	Toxic
50-00-0	Ignitable
Methyl alcohol	Toxic
67-56-1	Ignitable

14. TRANSPORT INFORMATION

Transportation information is provided as a general reference only and may not be applicable in all situations. This information applies to non-bulk shipments only. Per 49 CFR §173.22, it is the shipper's responsibility to ensure that all materials are properly packaged, classified and labeled prior to shipment.

DOT	Not regulated
	Not regulated

IATA Not regulated

15. REGULATORY INFORMATION

US Federal Regulations

<u>SARA 313</u>

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product contains a chemical or chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

Chemical Name	SARA 313 - Threshold Values %	
Formaldehyde - 50-00-0	0.1	

Methyl alcohol - 67-56-1	1.0	
SARA 311/312 Hazard Categories		
Acute health hazard	Yes	
Chronic Health Hazard	Yes	
Fire hazard	No	
Sudden release of pressure hazard	No	
Reactive Hazard	No	

CWA (Clean Water Act)

This product contains the following substances which are regulated pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
Formaldehyde 50-00-0	100 lb	-	-	х
Sodium Phosphate, dibasic 7558-79-4	5000 lb	-	-	Х

CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302)

Chemical Name	Hazardous Substances RQs	CERCLA/SARA RQ	Reportable Quantity (RQ)
Formaldehyde 50-00-0	100 lb	100 lb	RQ 100 lb final RQ RQ 45.4 kg final RQ
Methyl alcohol 67-56-1	5000 lb	-	RQ 5000 lb final RQ RQ 2270 kg final RQ
Sodium Phosphate, dibasic 7558-79-4	5000 lb	-	RQ 5000 lb final RQ RQ 2270 kg final RQ

US State Regulations

California Proposition 65

This product contains the following Proposition 65 chemicals

Chemical Name	California Proposition 65	
Formaldehyde - 50-00-0	Carcinogen	
Methyl alcohol - 67-56-1	Developmental	

U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Formaldehyde 50-00-0	Х	X	Х
Methyl alcohol 67-56-1	Х	X	Х
Sodium Phosphate, dibasic 7558-79-4	Х	X	Х

16. OTHER INFORMATION

Prepared By
Revision Date
Disclaimer

CDI Regulatory Affairs (Email: compliance@cancerdiagnostics.com) 2021-08-13

Cancer Diagnostics, Inc. makes no representation as to the comprehensiveness or accuracy of this document. Individuals using this information, or the product to which it refers, must exercise their independent judgment in determining all appropriateness for a particular purpose. Accordingly, Cancer Diagnostics, Inc. will not be responsible for damages of any kind or nature resulting from the use of this information or the corresponding product. No representations or warranties of any kind or nature, including but not limited to: express warranties, implied warranties or merchantability, or warranties of fitness for a particular purpose, are made hereunder with respect to the information set forth herein or to the product to which the information refers.

End of Safety Data Sheet



Addendum to Formalin NBF SDS

Item Codes for 10% Formalin Fixative, Neutral, Phosphate Buffered

FX0030, FX0055, FX1000, FX1001, FX1002, FX1003, FX1004, FX1004-GR, FX1005, FX1006, FX1006-576, FX1006- CSTM-1735, FX1006-CSTM-2834, FX1006-CSTM-301, FX1006-H, FX1006-O, FX1006-KT, FX1007, FX1007-KT, FX1007-H, FX1007- 128, FX1008, FX1008-H, FX1008-CSTM-2834, FX1008-CSTM- 301, FX1009, FX1009-H, FX1010, FX1010-CSTM-2834, FX1010- H, FX1010-KT, FX1011, FX1012, FX1013, FX0030, FX0055, 004096, 004096-200, PP1006, PP1008, PP1010, PP1012, PP1016, PP1012-25, PP1010-301, PP1012-301, PP1012-301-BX, PP1016-301, PP1016-BOX, PP1012-4438, FX1006-HPY, FX1007-HPY, FX1008-HPY, FX1009-HPY, FX1010-HPY, FX1011-24, FX1011-CUSTOM FILL-180-190, FX1116-24, FX1132-12, FX1006-BHZD, FX1007-BHZD, FX1008-BHZD, FX1009-BHZD, FX1010-BHZD, FX1000-20, FX1006-20FMA, FX1006-384, FX1006-O, FX1008-245, FX1006-CS, FX1007-CS, FX1008-CS, FX1009-CS, FX1010-CS, FX1004-GR, FX1006-100, FX1006-128, FX1006-BX, FX1007-BX, FX1008-BX, FX1009-BX, FX1010-BX, FX1006-HPY-768, FX1006-CSTM-1735, FX1006-KT-PPL, FX1007-CSTM-1735, FX1008-CSTM-2834, FX1008-CSTM-301

*Also Included, any item code listed above that may also be prefixed with "SAMPLE-"

1 | Page 2021-08-13

Tel: (877) 846-5393 | Fax: (877) 817-1716 | www.cancerdiagnostics.com