

1. Identification
Product Name: All-Fix™

Item #: FX1068, FX1069, FX1070, FX1071, FX1072, FX1073, FX1074, FX1075

Synonyms: N/A

Recommended Use: N/A

Manufacturer/Supplier:

 Cancer Diagnostics, Inc.
 116 Page Point Circle
 Durham, NC 27703
 1-877-846-5393

Restrictions on Use: N/A

In Case of Emergency:

 Chemtrec US 1-800-424-9300
 Infotrac International 1-352-323-3500

2. Hazards Identification
OSHA Hazard Classification(s):

 Acute Toxicity - Inhalation - Category 4
 Skin Irritation - Category 2
 Eye Irritation - Category 2A
 Specific Target Organ Toxicity (single exposure) - Category 1
 Specific Target Organ Toxicity (repeated exposure) - Category 2
 Toxic to Reproduction - Category 2
 Germ Cell Mutagenicity - Category 2
 Sensitization - Skin - Category 1B

Signal Word: Danger

Hazard Statement(s): Harmful if inhaled. Causes skin irritation. Causes serious eye irritation. Causes damage to organs(respiratory system). May cause damage to organs(respiratory system, central nervous system, liver, blood) through prolonged or repeated exposure. Suspected of damaging fertility of the unborn child. Suspect of causing genetic defects. May cause an allergic skin reaction.

Pictogram(s):

Precautionary Statement(s): Prevention: Avoid breathing dust, vapors. Use only outdoors or in a well-ventilated area. Wash body thoroughly after handling. Wear protective gloves. Wear eye protection, face protection. Do not breathe dust, vapors. Do not eat, drink or smoke when using this product. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves, protective clothing, eye protection and face protection. Contaminated work clothing must not be allowed out of the workplace.

Response: If inhaled: Remove person to fresh air and keep comfortable for breathing. Call a doctor if you feel unwell. If on skin: Wash with plenty of water. Specific treatment (see first aid section on this label). If skin irritation or rash occurs: Get medical attention. Take off all contaminated clothing and wash it before reuse. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing If eye irritation persists: Get medical attention. If exposed or concerned: Call a doctor. If exposed or concerned: Get medical attention.

Storage: Store locked up.

Disposal: Dispose of contents/container in accordance with local regulations.

Descriptions of Hazards not otherwise classified: N/A

Percent of mixture with unknown acute toxicity: N/A

3. Composition and Information on Ingredients

| Chemical Name | Common Name | CAS # | Concentration % |
|---------------|-------------|-----------|-----------------|
| Water | | 7732-18-5 | Trade Secret |
| Ethyl Alcohol | | 64-17-5 | Trade Secret |

| | | | |
|---------|--|----------|--------------|
| Glyoxal | | 107-22-2 | Trade Secret |
|---------|--|----------|--------------|

4. First Aid Measures

Eye Contact: If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Skin Contact: If on skin (or hair): Take off immediately all contaminated clothing and wash before reuse. Wash with plenty of water. If skin irritation occurs: Get medical advice/attention.

Inhalation: Remove to fresh air; give artificial respiration if breathing has stopped. Get medical advice/attention if you feel unwell.

Ingestion: Do not induce vomiting. If vomiting occurs spontaneously, keep head below hips to prevent aspiration of liquid into the lungs. Get medical attention.

Symptoms: Irritation eyes, nose, throat; headache, dizziness

Recommendations for immediate medical care/special treatment: Get medical advice/attention if you feel unwell.

5. Fire- Fighting Measures

Extinguishing Media: Dry chemical, carbon dioxide, alcohol foam, water.

Fire Hazards (Chemical): Not flammable.

Special Protective Equipment: Fire fighters should use self-contained breathing apparatus and protective clothing.

Precautions for Firefighters: Fire fighters should use self-contained breathing apparatus and protective clothing.

6. Accidental Release Measures

Emergency Procedures: Evacuate the area of all unnecessary personnel. Wear suitable protective equipment. Eliminate all sources of ignition and provide ventilation.

Protective Equipment: See section 8

Environmental Precautions: Prevent release to the environment by using barriers.

Containment and Clean-Up Procedures: Use barriers to prevent spreading. Collect spill in container. Call waste authorities.

7. Handling and Storage

Handling: Do not breathe vapors. Do not eat, drink or smoke when using this product.

Storage: Store locked up. Store in a well-ventilated, cool place. Keep container lid closed.

8. Exposure Controls/Personal Protection

OSHA Permissible Exposure Limits (PELs):

| Reagent | CAS # | OSHA PEL TWA |
|---------------|---------|--------------|
| Ethyl Alcohol | 64-17-5 | 1000ppm |

ACGIH Threshold Limit Values (TLVs):

| Reagent | CAS # | ACGIH PEL TLV | ACGIH STEL |
|---------------|---------|---------------|------------|
| Ethyl Alcohol | 64-17-5 | 1000ppm | 1000ppm |

Engineering Controls: Use in a well ventilated area to prevent exposure. Maintain eyewash fountain and quick-drench facilities in work areas.

Personal Protective Measures: Wear gloves, lab coat, eye protection and impervious footwear. Contact lenses should not be worn when working with this material.

Special PPE Requirements: If ventilation hood not available wear respirator.

9. Physical and Chemical Properties Section

Appearance: Colorless, Liquid

Molecular Weight: N/A
Molecular Formula: N/A
pH: 5.6-5.8
Boiling Point and Boiling Range: N/A
Melting Point/Freezing Point: N/A
Flash Point: N/A
Specific Gravity/Relative Density: N/A
Odor: N/A
Odor Threshold: N/A
Color: Colorless
Flammability (solid/gas): N/A
Vapor Density: N/A
Upper/Lower flammability or explosive limits: N/A
Vapor Pressure: N/A
Evaporation Rate: N/A
Partition Coefficient: n-octanol/water: N/A
Viscosity: N/A
Auto-ignition temperature: N/A
Solubility: Miscible in water
Decomposition Temperature: N/A

10. Stability and Reactivity

Reactivity: N/A
Chemical Stability: Stable
Conditions of Stability/Instability: Stable
Stabilizers needed: None
Safety issue indicated by appearance change: N/A
Other: N/A
Hazardous Reactions: N/A
Hazardous Polymerization: Does not occur
Conditions to avoid: N/A
Classes of Incompatible Materials: Oxidizers, Strong Acids, Strong Bases
Hazardous Decomposition Products: Thermal-oxidation degradation can produce oxides of carbon. Toxic gases and vapors (I.e. Carbon monoxide) may be released in a fire.

11. Toxicological Information

Likely Routes of Exposure

Eyes: Irritation. May cause permanent damage.
Skin: Irritation. May cause chemical burns or dermatitis.
Inhalation: Vapors are harmful by inhalation. Dizziness, headache, nausea. Irritation to respiratory tract or lungs.
Ingestion: Nausea.

Signs or Symptoms of Exposure: Nausea.

Effects from short term exposure (delayed, immediate, chronic): Irritation to the eyes, nose, throat; headache, dizziness, nausea.

Acute Toxicity (Numerical Measures): N/A

Carcinogenicity (NTP, IARC, OSHA): N/A

12. Ecological Information

Ecotoxicity:
Persistence and degradability:
Bioaccumulation Potential (octanol-water partition coefficient, BCF):
Mobility in the soil:
Adverse Environmental Effects:

13. Disposal Considerations

Recommended Disposal Containers: Check with your local waste authorities*
Recommended Disposal Methods: Do not dispose of in drains, check with your local waste authorities.*
Physical/Chemical Properties affecting Disposal: See section 2 and section 9 applicable information.*
Special Precautions for Landfill and Incineration Activities: Check with your local waste authorities.*
Waste Stream: Consult your local or regional authorities.*

14. Transport Information

| | |
|---|--|
| 14.1 UN Number DOT, IATA,IMDG, ADR | DOT: Not Regulated for Transport by DOT. ADR/RID: Not Regulated for Transport by ADR/RID. IMDG: Not Regulated for Transport by IMDG. IATA: Not Regulated for Transport by IATA. |
| 14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR | ADR/RID: Not Regulated for Transport by ADR/RID. IMDG: Not Regulated for Transport by IMDG. IATA: Not Regulated for Transport by IATA. DOT: Not Regulated for Transport by DOT. |
| 14.3 Transport Hazard Class(es) | DOT: Not Regulated for Transport by DOT. IATA: Not Regulated for Transport by IATA. IMDG: Not Regulated for Transport by IMDG. ADR/RID: Not Regulated for Transport by ADR/RID. |
| 14.4 Packing Group DOT, IATA,IMDG, ADR | ADR/RID: Not Regulated for Transport by ADR/RID. IMDG: Not Regulated for Transport by IMDG. IATA: Not Regulated for Transport by IATA. DOT: Not Regulated for Transport by DOT. |
| 14.5 Environmental Hazards | Marine Pollutant: No |
| 14.6 Special Precautions for User | Not applicable. |

15. Regulatory Information

OSHA:
DOT:
EPA:
CPSC:

16. Other Information


Revision Date: 2019-01-11

NFPA

| | |
|-----------------|---|
| Health | 2 |
| Fire Hazard | 0 |
| Reactivity | 0 |
| Specific Hazard | |

National Fire Protection Association (USA) NFPA

Fire Hazard



HMIS

| | |
|---------------------|---|
| Health | 2 |
| Flammability | 0 |
| Physical Hazard | 0 |
| Personal Protection | |

Hazardous Material Information System HMIS

| | |
|---------------------|---|
| Health | 2 |
| Flammability | 0 |
| Physical Hazard | 0 |
| Personal Protection | |

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.