

SAFETY DATA SHEET

Cancer Diagnostics, Inc.
 116 Page Point Circle
 Durham, NC 27703

Emergency Phone: 800-424-9300 4300
 Information Phone: 877-846-5393

SECTION I PRODUCT IDENTIFICATION

STYLE: VP0810, VP01114, VP1620, VP06014

SECTION II HAZARD IDENTIFICATION

None

SECTION III INGREDIENTS

	<u>CAS #:</u>	<u>OSHA PEL</u>	<u>ACGIN TLV</u>
2-Polyester	25038-59-9	NA	NA
3-Delnet P530 (Polyethylene)		NA	NA
4-Proprietary coupling agent		NA	NA

SECTION IV PHYSICAL & CHEMICAL PROPERTIES

PHYSICAL FORM	- Needlepunched Nonwoven
COLOR	-White / Blue Top Surface
ODOR	- N/A
PH	- Neutral
BOILING POINT	- N/A
MELTING POINT	- 1) NA 2) 485F 3) 275F
SOLUBILITY IN WATER	- N/A
BULK DENSITY	-10.3 oz/yd ²
VAPOR PRESSURE	- N/A

SECTION V

FIRE AND EXPLOSION DATA

FLASH POINT	- N/A
EXTINGUISHING MEDIA	- Water, Dry Chemical, Foam, Carbon Dioxide (CO ²)
FLAMMABILITY LIMITS	- N/A
AUTOIGNITION	- None
OTHER	- Material will burn when exposed to a flame

SECTION VI

FIRE FIGHTING MEASURES

SPECIAL FIRE FIGHTING PROCEDURES	- Keep personnel removed and upwind of fire. Wear self-contained breathing apparatus. Wear full protective equipment (full Bunker gear).
----------------------------------	--

SECTION VII

TOXICOLOGICAL INFORMATION

STABILITY	- There are no known health hazards associated with the use of this product.
SIGNS AND SYMPTOMS OF EXPOSURE	- None Known
CARCINOGENICITY	- The components of this fabric are not listed by NTP, IRAC or regulated as a carcinogen by OSHA
MEDICAL CONDITIONS AGGRAVATED BY EXPOSURE	- None expected or known.

SECTION VIII

STABILITY & REACTIVITY DATA

STABILITY	- Stable
POLYMERIZATION	- Will not occur
INCOMPATIBILITIES	- Strong acids, alkali, & oxidizers
INSTABILITY CONDITIONS DECOMPOSITION PRODUCTS	- NA

SECTION IX	HANDLING & STORAGE
-------------------	-------------------------------

SECTION X	ACCIDENTIAL RELEASE MEASURES
------------------	-------------------------------------

If the fabric is not contaminated, place into clean container for use. If contaminated, place into containers for proper disposal.

SECTION XI	DISPOSAL CONSIDERATIONS
-------------------	--------------------------------

Dispose of fabric in a licensed landfill facility according to applicable state and local laws, federal industrial waste regulations and other applicable standards. Incineration in a licensed facility is not recommended.

OTHER PRECAUTIONS

Keep work and storage areas clean
 Keep product away from flame

SECTION XII	EXPOSURE CONTROLS & PERSONAL PROTECTION
--------------------	--

EYE PROTECTION	None
SKIN PROTECTION	Good personal hygiene measures should include personal protective equipment to reduce exposure when extensive handling of fabric is necessary
RESPIRATORY REQUIREMENTS	N/A
VENTILATION REQUIREMENTS	Local exhaust is recommended
WORK / HYGIENIC PRACTICES	Wash hands and face before eating, drinking, smoking and using restroom facilities.

SECTION XIII	FIRST-AID MEASURES
---------------------	---------------------------

EYE	If fiber gets into eyes rinse with cold water if Irritation continues seek medical advice
SKIN	If irritation occurs rinse with cold water and use protective clothing; if irritation continues seek medical advice.

SECTION XIV
TRANSPORTATION

There are no restrictions for transportation; product does not contain any hazard substances.

SECTION XV
ECOLOGICAL INFORMATION

Material is not biodegradable; dispose of in accordance with local regulations

SECTION XVI
REGULATORY INFORMATION

None

SECTION
OTHER INFORMATION

REASON FOR ISSUE	Update
PREPARED BY	Compliance@cancerdiagnostics.com
APPROVED DATE	September 24, 2015
SUPERSEDES DATE	N/A
MSDS NUMBER	N/A

This information is furnished without warranty, expressed or implied, except that it is accurate to the best knowledge of Cancer Diagnostics, Inc.

Cancer Diagnostics assumes no legal responsibility for use or reliance upon this information.

HMIS RATING:	Health	Fire	Reactivity	Personal Protective Equipment (B)
Health	0	1	0	