



SAFETY DATA SHEET

Revision Date 2021-06-21

Version 1

1. IDENTIFICATION

Product Name PureView® Focus-AQ™ Concentrate

Product Code 16703-C, 16703-C4

Recommended Use For laboratory, scientific, R&D or manufacturing use.

Company Cancer Diagnostics, Inc.
116 Page Point Circle
Durham, NC 27703
Tel. (877) 846-5393

Emergency Telephone Call CHEMTREC 1-800-424-9300

2. HAZARDS IDENTIFICATION

Classification

OSHA Regulatory Status

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Label elements

Although this product does not normally require specific hazard precautions, chemical users should always take care to minimize personnel exposure and workplace contamination. Store tightly closed at 15-30C.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS No.	Weight-%
Water	7732-18-5	>98
Magnesium sulfate, heptahydrate	10034-99-8	<25.5
DL-Malic Acid	6915-15-7	<9
Sodium Citrate, Dihydrate	6132-04-3	<2
Sodium hydroxide	1310-73-2	>0.4

Chemical Name	CAS No.	Weight-%
Water	7732-18-5	>98
7732-18-5		
Magnesium sulfate, heptahydrate	10034-99-8	<25.5
10034-99-8		
DL-Malic Acid	6915-15-7	<9
6915-15-7		
Sodium Citrate, Dihydrate	6132-04-3	<2
6132-04-3		
Sodium hydroxide	1310-73-2	>0.4
1310-73-2		

4. FIRST AID MEASURES

Description of first aid measures

General advice This chemical is not expected to produce any significant adverse health effects.

Eye contact	Flush eyes with plenty of water, removing contact lenses if present. Get medical attention if irritation develops.
Skin contact	Wash thoroughly with soap and water while removing contaminated garments. Get medical attention if irritation develops.
Inhalation	Remove to fresh air. Get medical attention for any breathing difficulty.
Ingestion	Rinse mouth and drink several glasses of water. Contact a physician or poison control center if symptoms develop.

5. FIRE-FIGHTING MEASURES

Suitable extinguishing media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment

Specific hazards arising from the chemical

No Data Available.

Protective equipment and precautions for firefighters

Firefighters should wear self-contained breathing apparatus with full facepiece operated in pressure-demand or other positive pressure mode.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal precautions Use personal protective equipment as required.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Pick up and transfer to properly labeled containers.

7. HANDLING AND STORAGE

Precautions for safe handling Handle in accordance with good industrial hygiene and safety practice.

Storage Conditions Keep containers tightly closed in a dry, cool and well-ventilated place.

Incompatible materials None known based on information supplied.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Occupational exposure limits Not applicable
Appropriate engineering controls

Engineering Controls Emergency showers, eyewash stations, ventilation systems.

General Hygiene Considerations Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical state Liquid
Appearance Clear liquid

Odor	No Data Available
Odor threshold	No Data Available
pH	<5
Melting point / freezing point	No Data Available
Boiling point / boiling range	No Data Available
Flash point	No Data Available
Evaporation rate	No Data Available
Flammability (solid, gas)	No Data Available
Flammability Limit in Air	
Upper flammability limit:	No Data Available
Lower flammability limit:	No Data Available
Vapor pressure	No Data Available
Vapor density	No Data Available
Relative density	No Data Available
Water solubility	Miscible with water
Solubility in other solvents	No Data Available
Partition coefficient	No Data Available
Autoignition temperature	No Data Available
Decomposition temperature	No Data Available
Kinematic viscosity	No Data Available

10. STABILITY AND REACTIVITY

Stability	Stable under recommended storage conditions.
Possibility of Hazardous Reactions	No Data Available.
Conditions to avoid	Extremes of temperature and direct sunlight
Incompatible materials	None known based on information supplied.
Hazardous Decomposition Products	None known based on information supplied.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Component	Oral LD50	Dermal LD50	Inhalation LC50
-	> 90 mL/kg (Rat)	-	-
-	= 325 mg/kg (Rat)	= 1350 mg/kg (Rabbit)	-

Symptoms	No Data Available.
----------	--------------------

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation	No Data Available.
Serious eye damage/eye irritation	No Data Available.
Irritation	No Data Available.
Sensitization	No Data Available.
Germ cell mutagenicity	No Data Available.
Carcinogenicity	This chemical does not contain any carcinogens or potential carcinogens as listed by ACGIH, OSHA, IARC or NTP

12. ECOLOGICAL INFORMATION

Ecotoxicity

No information available

Persistence and degradability

No Data Available.

Bioaccumulation

No Data Available.

13. DISPOSAL CONSIDERATIONS**Waste treatment methods****Disposal of wastes**

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Chemical Name	California Hazardous Waste Status
Sodium hydroxide 1310-73-2	Toxic Corrosive

14. TRANSPORT INFORMATION

Note: Transportation information is provided as a general reference only and may not be applicable in all situations. This information applies to non-bulk shipments only. Per 49 CFR §173.22, it is the shipper's responsibility to ensure that all materials are properly packaged, classified and labeled prior to shipment.

DOT

Not regulated

IATA

Not regulated

15. REGULATORY INFORMATION

No Data Available.

16. OTHER INFORMATION**Prepared By**CDI Regulatory Affairs (Email: compliance@cancerdiagnostics.com)**Revision Date**

2021-06-21

Disclaimer

Cancer Diagnostics, Inc. makes no representation as to the comprehensiveness or accuracy of this document. Individuals using this information, or the product to which it refers, must exercise their independent judgment in determining all appropriateness for a particular purpose. Accordingly, Cancer Diagnostics, Inc. will not be responsible for damages of any kind or nature resulting from the use of this information or the corresponding product. No representations or warranties of any kind or nature, including but not limited to: express warranties, implied warranties or merchantability, or warranties of fitness for a particular purpose, are made hereunder with respect to the information set forth herein or to the product to which the information refers.

End of Safety Data Sheet