

Hydrochloric Acid, 20%

1. Identification

Product Name: Hydrochloric Acid, 20%

Item #: SSC1085

Synonyms: N/A

Recommended Use: Special Stains, Laboratory Reagent

Restrictions on Use: Any use other than recommended.

Manufacturer:

Cancer Diagnostics, Inc.

116 Page Point Circle

Durham, NC 27703

1-877-846-5393

In Case of Emergency:

Chemtrec US 1-800-424-9300

Infotrac International 1-352-323-3500

2. Hazards Identification

OSHA Hazard Classification(s):

Acute Toxicity - Inhalation - Category 4

Skin Corrosion - Category 1B

Eye Damage - Category 1

Signal Word: Danger

Hazard Statement(s): Harmful if inhaled. Causes severe skin burns and eye damage. Causes serious eye damage.

Pictogram(s):

Precautionary Statement(s): Prevention: Avoid breathing dust, vapors. Use only outdoors or in a well-ventilated area. Do not breathe dusts or mists. Wash body thoroughly after handling. Wear protective gloves, protective clothing, eye protection and face protection. Wear eye protection, face protection.

Response: If inhaled: Remove person to fresh air and keep comfortable for breathing. Call a doctor if you feel unwell. If swallowed: Rinse mouth. Do NOT induce vomiting. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. Take off all contaminated clothing and wash it before reuse. Immediately call a doctor. Specific treatment (see first aid section on this label). If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing

Storage: Store locked up.

Disposal: Dispose of contents/container in accordance with local regulations.

Descriptions of Hazards not otherwise classified: N/A

Percent of mixture with unknown acute toxicity: N/A

3. Composition and Information on Ingredients

Chemical Name	Common Name	CAS #	Concentration %
Hydrochloric Acid	HCl	7647-01-0	20
Water		7732-18-5	80

4. First Aid Measures

Eye Contact: If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Skin Contact: If on skin (or hair): Take off immediately all contaminated clothing and wash before reuse. Wash with plenty of water. If skin irritation occurs: Get medical advice/attention.

Inhalation: Remove to fresh air; give artificial respiration if breathing has stopped. Get medical advice/attention if you feel unwell.

Ingestion: Do not induce vomiting. If vomiting occurs spontaneously, keep head below hips to prevent aspiration of liquid into the lungs. Get medical attention.

Hydrochloric Acid, 20%

Symptoms: Irritation eyes, nose, throat; headache, dizziness. Irritation to skin.

Recommendations for immediate medical care/special treatment: Get medical advice/attention if you feel unwell.

5. Fire- Fighting Measures

Extinguishing Media: Dry chemical, carbon dioxide, alcohol foam, water.

Fire Hazards (Chemical): Not flammable.

Special Protective Equipment: Fire fighters should use self-contained breathing apparatus and protective clothing.

Precautions for Firefighters: Carbon monoxide and unidentified organic compounds may be formed during combustion.

6. Accidental Release Measures

Emergency Procedures: Evacuate the area of all unnecessary personnel. Wear suitable protective equipment. Eliminate all sources of ignition and provide ventilation.

Protective Equipment: See section 8

Environmental Precautions: Prevent release to the environment by using barriers.

Containment and Clean-Up Procedures: Use barriers to prevent spreading. Collect spill in container. Call waste authorities.

7. Handling and Storage

Handling: Do not breathe vapors. Do not eat, drink or smoke when using this product.

Storage: Store locked up. Store in a cool, well-ventilated place. Keep lid tightly sealed when not in use.

8. Exposure Controls/Personal Protection

OSHA Permissible Exposure Limits (PELs):

Reagent	CAS #	OSHA PEL TWA
Hydrochloric Acid	7647-01-0	5 ppm (7 mg/m ³) Ceiling

ACGIH Threshold Limit Values (TLVs):

Reagent	CAS #	ACGIH PEL TLV	ACGIH STEL
Hydrochloric Acid	7647-01-0	2 ppm (3 mg/m ³) Ceiling	

Engineering Controls: Use in a well ventilated area to prevent exposure. Maintain eyewash fountain and quick-drench facilities in work areas.

Personal Protective Measures: Wear gloves, lab coat, eye protection and impervious footwear. Contact lenses should not be worn when working with this material.

Special PPE Requirements: If ventilation hood not available wear respirator.

9. Physical and Chemical Properties Section

Appearance: Colorless, Liquid

Molecular Weight: N/A

Molecular Formula: N/A

pH: Less than 2

Boiling Point and Boiling Range: N/A

Melting Point/Freezing Point: N/A

Flash Point: N/A

Specific Gravity/Relative Density: N/A

Odor: N/A

Odor Threshold: N/A

Color: Colorless

Flammability (solid/gas): N/A

Hydrochloric Acid, 20%

Vapor Density: N/A
Upper/Lower flammability or explosive limits: N/A
Vapor Pressure: N/A
Evaporation Rate: N/A
Partition Coefficient: n-octanol/water: N/A
Viscosity: N/A
Auto-ignition temperature: N/A
Solubility: N/A
Decomposition Temperature: N/A

10. Stability and Reactivity

Reactivity:
Chemical Stability: Stable
Conditions of Stability/Instability: Stable under normal conditions of temperature and pressure.
Stabilizers needed: None
Safety issue indicated by appearance change: N/A
Other: N/A
Hazardous Reactions: N/A
Hazardous Polymerization: Does not occur
Conditions to avoid: N/A
Classes of Incompatible Materials: Oxidizers, Strong Acids, Strong Bases
Hazardous Decomposition Products: Thermal-oxidation degradation can produce oxides of carbon. Toxic gases and vapors (i.e. Carbon monoxide) may be released in a fire.

11. Toxicological Information

Likely Routes of Exposure

Eyes: Corrosive to eyes. May cause irritation, reddening and/or permanent damage. May cause loss of vision.
Skin: Corrosive to skin. May cause irritation, reddening, dermatitis and/or permanent scars.
Inhalation: Dizziness, headache. May cause respiratory distress if a large amount is inhaled. Vapors are strong and cause coughing spells when inhaled. May trigger asthma attacks or asthma like symptoms in individuals with respiratory sensitivity.
Ingestion: Corrosive solution will cause chemical burns to mouth, throat and stomach. Ingestion of large amounts may be fatal.

Signs or Symptoms of Exposure: Nausea, excessive coughing, asthma like symptoms, headaches, dermatitis or reddening of skin.

Effects from short term exposure (delayed, immediate, chronic): Irritation to the eyes, nose, throat; headache, dizziness, nausea. Irritation to respiratory system.

Acute Toxicity (Numerical Measures): Concentrated Hydrochloric Acid: LD50(oral, rat)=900 mg/kg; LC50(inhalation, mouse)=1108 ppm/1H; LC50(inhalation, mouse)=3940 mg/m³/30M

Carcinogenicity (NTP, IARC, OSHA): Not listed as a carcinogen.

12. Ecological Information

Ecotoxicity: Concentrated Hydrochloric Acid CAS 7647-01-0

Fish: LC50 (96 Hr) Mosquito Fish: 282 mg/L, LC100(24Hr) Trout: 10mg/L

Invertebrates: LC50(48Hr) Starfish: 100-330 mg/L, LC50 (48Hr)

Shrimp: 100-330 mg/L

Persistence and degradability: N/A

Bioaccumulation Potential (octanol-water partition coefficient, BCF): N/A

Mobility in the soil: N/A

Hydrochloric Acid, 20%

Adverse Environmental Effects: N/A

13. Disposal Considerations

Recommended Disposal Containers: Check with your local waste authorities.*

Recommended Disposal Methods: Do not dispose of in drains, check with your local waste authorities.*

Physical/Chemical Properties affecting Disposal: See section 2 and section 9 applicable information.*

Special Precautions for Landfill and Incineration Activities: Check with your local waste authorities.*

Waste Stream: Consult your local or regional authorities.*

14. Transport Information

UN Number: UN1789

UN Proper Shipping Name: Hydrochloric Acid

Transport Hazard Class(es): 8

Packing Group Number: II

Environmental Hazards (IMDG code):

Marine Pollutant: No

Transport in Bulk (IBC Code): N/A

Special Transport Precautions: N/A

15. Regulatory Information

OSHA:

DOT:

EPA:

CPSC:

Hydrochloric Acid, 20%

16. Other Information


Revision Date: 2020-11-13

NFPA

Health	3
Fire Hazard	0
Reactivity	2
Specific Hazard	COR

National Fire Protection Association (USA) NFPA

Fire Hazard



Specific Hazard

HMIS

Health	3
Flammability	0
Physical Hazard	2
Personal Protection	

Hazardous Material Information System HMIS

Health	3
Flammability	0
Physical Hazard	2
Personal Protection	

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Revises: Revision ID:
2019-01-11 SDS1003