

OG-6 Powdered Stain

1. Identification

Product Name: OG-6 Powdered Stain

Item #: CMG001, CMG025, CMG050

Synonyms: N/A

Recommended Use: Laboratory reagent. Stain

Restrictions on Use: Any use other than recommended.

Manufacturer:

Cancer Diagnostics, Inc.
116 Page Point Circle
Durham, NC 27703
1-877-846-5393

In Case of Emergency:

Chemtec US 1-800-424-9300
Infotrac International 1-352-323-3500

2. Hazards Identification

OSHA Hazard Classification(s):

No OSHA Hazard Classifications Applicable

Signal Word: N/A

Hazard Statement(s): N/A

Pictogram(s): N/A

Precautionary Statement(s): Prevention: N/A

Response: N/A

Storage: N/A

Disposal: N/A

Descriptions of Hazards not otherwise classified: N/A

Percent of mixture with unknown acute toxicity: N/A

3. Composition and Information on Ingredients

Chemical Name	Common Name	CAS #	Concentration %
OG-6 Stain(s)		Trade Secret	Trade Secret

4. First Aid Measures

Eye Contact: Flush eyes with water as a precaution.

Skin Contact: Wash off with plenty of water. Remove contaminated clothing and launder before reuse as a precaution.

Inhalation: Move person to fresh air; give artificial respiration if breathing has stopped.

Ingestion: Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if discomfort occurs.

Symptoms: N/A

Recommendations for immediate medical care/special treatment: If exposure by any route causes irritation get medical advice/attention.

5. Fire- Fighting Measures

Extinguishing Media: Dry chemical, carbon dioxide, alcohol foam, water.

Fire Hazards (Chemical): Not flammable.

Special Protective Equipment: Fire fighters should use self-contained breathing apparatus and protective clothing.

Precautions for Firefighters: Carbon monoxide and unidentified organic compounds may be formed during combustion.

6. Accidental Release Measures

Emergency Procedures: Evacuate the area of all unnecessary personnel. Wear suitable protective equipment. Eliminate all sources of ignition and provide ventilation.

OG-6 Powdered Stain

Protective Equipment: See section 8

Environmental Precautions: Prevent release to the environment by using barriers.

Containment and Clean-Up Procedures: Use barriers to prevent spreading. Collect spill in container. Call waste authorities.

7. Handling and Storage

Handling: Do not breathe vapors. Do not eat, drink or smoke when using this product.

Storage: Store in a well-ventilated place. Keep cool. Keep lid tightly closed.

8. Exposure Controls/Personal Protection

OSHA Permissible Exposure Limits (PELs):

Reagent	CAS #	OSHA PEL TWA
N/A		

ACGIH Threshold Limit Values (TLVs):

Reagent	CAS #	ACGIH PEL TLV	ACGIH STEL
N/A			

Engineering Controls: Use in a well ventilated area to prevent exposure. Maintain eyewash fountain and quick-drench facilities in work areas.

Personal Protective Measures: Wear gloves, lab coat, eye protection and impervious footwear. Contact lenses should not be worn when working with this material.

Special PPE Requirements: If ventilation hood not available wear respirator.

9. Physical and Chemical Properties Section

Appearance: Orange, Powder

Molecular Weight: N/A

Molecular Formula: N/A

pH: N/A

Boiling Point and Boiling Range: N/A

Melting Point/Freezing Point: N/A

Flash Point: N/A

Specific Gravity/Relative Density: N/A

Odor: N/A

Odor Threshold: N/A

Color: Orange

Flammability (solid/gas): N/A

Vapor Density: N/A

Upper/Lower flammability or explosive limits: N/A

Vapor Pressure: N/A

Evaporation Rate: N/A

Partition Coefficient: n-octanol/water: N/A

Viscosity: N/A

Auto-ignition temperature: N/A

Solubility: N/A

Decomposition Temperature: N/A

10. Stability and Reactivity

Reactivity:

OG-6 Powdered Stain

Chemical Stability: Stable

Conditions of Stability/Instability: Stable under normal conditions of temperature and pressure.

Stabilizers needed: None

Safety issue indicated by appearance change: N/A

Other: N/A

Hazardous Reactions: N/A

Hazardous Polymerization: Does not occur

Conditions to avoid: Incompatible materials, dust generation, excess heat.

Classes of Incompatible Materials: Oxidizers, Strong Acids, Strong Bases

Hazardous Decomposition Products: Hydrogen chloride, carbon monoxide, oxides of nitrogen, oxides of sulfur, irritating and toxic fumes and gases, carbon dioxide

11. Toxicological Information

Likely Routes of Exposure

Eyes: Irritation.

Skin: Irritation.

Inhalation: Dizziness, headache.

Ingestion: Nausea.

Signs or Symptoms of Exposure: Nausea.

Effects from short term exposure (delayed, immediate, chronic): Irritation to the eyes, nose, throat; headache, dizziness, nausea.

Acute Toxicity (Numerical Measures): N/A

Carcinogenicity (NTP, IARC, OSHA): IARC: Group 3 carcinogen.

12. Ecological Information

Ecotoxicity: N/A

Persistence and degradability: N/A

Bioaccumulation Potential (octanol-water partition coefficient, BCF): N/A

Mobility in the soil: N/A

Adverse Environmental Effects: N/A

13. Disposal Considerations

Recommended Disposal Containers: Check with your local waste authorities*

Recommended Disposal Methods: Do not dispose of in drains, check with your local waste authorities.*

Physical/Chemical Properties affecting Disposal: See section 2 and section 9 applicable information.*

Special Precautions for Landfill and Incineration Activities: Check with your local waste authorities.*

Waste Stream: Consult your local or regional authorities.*

OG-6 Powdered Stain

14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not Regulated for Transport by DOT. ADR/RID: Not Regulated for Transport by ADR/RID. IMDG: Not Regulated for Transport by IMDG. IATA: Not Regulated for Transport by IATA.
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not Regulated for Transport by ADR/RID. IMDG: Not Regulated for Transport by IMDG. IATA: Not Regulated for Transport by IATA. DOT: Not Regulated for Transport by DOT.
14.3 Transport Hazard Class(es)	DOT: Not Regulated for Transport by DOT. IATA: Not Regulated for Transport by IATA. IMDG: Not Regulated for Transport by IMDG. ADR/RID: Not Regulated for Transport by ADR/RID.
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not Regulated for Transport by ADR/RID. IMDG: Not Regulated for Transport by IMDG. IATA: Not Regulated for Transport by IATA. DOT: Not Regulated for Transport by DOT.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

15. Regulatory Information

OSHA: N/A

DOT: N/A

EPA: N/A

CPSC: N/A

OG-6 Powdered Stain


16. Other Information

Revision Date: 2019-01-11

NFPA

Health	1
Fire Hazard	0
Reactivity	0
Specific Hazard	


National Fire Protection Association (USA) NFPA



HMIS

Health	1
Flammability	0
Physical Hazard	0
Personal Protection	

Hazardous Material Information System HMIS



Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.