

## Harris Hematoxylin

### 1. Identification

**Product Name:** Harris Hematoxylin

**Item #:** SH3777, SH4777, SH5777

**Synonyms:** Harris Hematoxylin without Mercuric Chloride

**Recommended Use:** N/A

**Restrictions on Use:** N/A

**Manufacturer:**
**In Case of Emergency:**

Cancer Diagnostics, Inc.

Chemtrec US 1-800-424-9300

116 Page Point Circle

Infotrac International 1-352-323-3500

Durham, NC 27703

1-877-846-5393

### 2. Hazards Identification

**OSHA Hazard Classification(s):**

Skin Irritation - Category 2

Eye Irritation - Category 2A

Specific Target Organ Toxicity (repeated exposure) - Category 2

Toxic to Reproduction - Category 2

**Signal Word:** Warning

**Hazard Statement(s):** Causes skin irritation. Causes serious eye irritation. May cause damage to organs(kidney) through prolonged or repeated exposure. Suspected of damaging fertility of the unborn child.

**Pictogram(s):**

**Precautionary Statement(s):** Prevention: Wash body thoroughly after handling. Wear protective gloves. Wear eye protection, face protection. Do not breathe dust, vapors. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves, protective clothing, eye protection and face protection.

Response: If on skin: Wash with plenty of water. Specific treatment (see first aid section on this label). If skin irritation or rash occurs: Get medical attention. Take off all contaminated clothing and wash it before reuse. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing If eye irritation persists: Get medical attention. Call a doctor if you feel unwell. If exposed or concerned: Get medical attention.

Storage: Store locked up.

Disposal: Dispose of contents/container in accordance with local regulations.

**Descriptions of Hazards not otherwise classified:** N/A

**Percent of mixture with unknown acute toxicity:** N/A

### 3. Composition and Information on Ingredients

Chemical Name	Common Name	CAS #	Concentration %
Water		732-18-5	Trade Secret
Hematoxylin		517-28-2	Trade Secret
Aluminum Salts		Trade Secret	Trade Secret
Select Glycols		Trade Secret	Trade Secret
Select Acid		Trade Secret	Trade Secret
Select Preservative		Trade Secret	Trade Secret

### 4. First Aid Measures

**Eye Contact:** If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

## Harris Hematoxylin

**Skin Contact:** If on skin (or hair): Take off immediately all contaminated clothing and wash before reuse. Wash with plenty of water.

If skin irritation occurs: Get medical advice/attention.

**Inhalation:** Remove to fresh air; give artificial respiration if breathing has stopped. Get medical advice/attention if you feel unwell.

**Ingestion:** Do not induce vomiting. If vomiting occurs spontaneously, keep head below hips to prevent aspiration of liquid into the lungs. Get medical attention.

**Symptoms:** Irritation eyes, nose, throat; headache, dizziness

**Recommendations for immediate medical care/special treatment:** Get medical advice/attention if you feel unwell.

### 5. Fire- Fighting Measures

**Extinguishing Media:** Dry chemical, carbon dioxide, alcohol foam, water.

**Fire Hazards (Chemical):** Not flammable.

**Special Protective Equipment:** Fire fighters should use self-contained breathing apparatus and protective clothing.

**Precautions for Firefighters:** Carbon monoxide and unidentified organic compounds may be formed during combustion.

### 6. Accidental Release Measures

**Emergency Procedures:** Evacuate the area of all unnecessary personnel. Wear suitable protective equipment. Eliminate all sources of ignition and provide ventilation.

**Protective Equipment:** See section 8

**Environmental Precautions:** Prevent release to the environment by using barriers.

**Containment and Clean-Up Procedures:** Use barriers to prevent spreading. Collect spill in container. Call waste authorities.

### 7. Handling and Storage

**Handling:** Do not breathe vapors. Do not eat, drink or smoke when using this product.

**Storage:** Store locked up. Keep container tightly closed. Store away from heat, sparks, open flame. Store in a cool, well-ventilated place.

### 8. Exposure Controls/Personal Protection

**OSHA Permissible Exposure Limits (PELs):**

Reagent	CAS #	OSHA PEL TWA
Glacial Acetic Acid	64-19-7	10 ppm, 25 mg/m3

**ACGIH Threshold Limit Values (TLVs):**

Reagent	CAS #	ACGIH PEL TLV	ACGIH STEL
Glacial Acetic Acid	64-19-7	10 ppm, 25 mg/m3	15 ppm, 37 mg/m3

**Engineering Controls:** Use in a well ventilated area to prevent exposure. Maintain eyewash fountain and quick-drench facilities in work areas.

**Personal Protective Measures:** Wear gloves, lab coat, eye protection and impervious footwear. Contact lenses should not be worn when working with this material.

**Special PPE Requirements:** If ventilation hood not available wear respirator.

### 9. Physical and Chemical Properties Section

**Appearance:** Dark Purple, Liquid

**Molecular Weight:** N/A

**Molecular Formula:** N/A

**pH:** 2.3-2.5

**Boiling Point and Boiling Range:** N/A

**Melting Point/Freezing Point:** N/A

**Flash Point:** N/A

## Harris Hematoxylin

**Specific Gravity/Relative Density:** N/A  
**Odor:** N/A  
**Odor Threshold:** N/A **Color:**  
Dark Purple **Flammability**  
**(solid/gas):** N/A **Vapor**  
**Density:** N/A  
**Upper/Lower flammability or explosive limits:** N/A  
**Vapor Pressure:** N/A  
**Evaporation Rate:** N/A  
**Partition Coefficient: n-octanol/water:** N/A  
**Viscosity:** N/A  
**Auto-ignition temperature:** N/A  
**Solubility:** N/A  
**Decomposition Temperature:** N/A

### 10. Stability and Reactivity

**Reactivity:**  
**Chemical Stability:** Stable  
**Conditions of Stability/Instability:** N/A  
**Stabilizers needed:** None  
**Safety issue indicated by appearance change:** N/A  
**Other:** N/A  
**Hazardous Reactions:** N/A  
**Hazardous Polymerization:** Does not occur  
**Conditions to avoid:** N/A  
**Classes of Incompatible Materials:** Oxidizers, Strong Acids, Strong Bases  
**Hazardous Decomposition Products:** Thermal-oxidation degradation can produce oxides of carbon. Toxic gases and vapors (I.e. Carbon monoxide) may be released in a fire.

### 11. Toxicological Information

#### Likely Routes of Exposure

**Eyes:** Irritation.  
**Skin:** Irritation.  
**Inhalation:** Dizziness, headache.  
**Ingestion:** Nausea. Dizziness.

**Signs or Symptoms of Exposure:** Nausea.

**Effects from short term exposure (delayed, immediate, chronic):** Irritation to the eyes, nose, throat; headache, dizziness, nausea.

**Acute Toxicity (Numerical Measures):** Glacial Acetic Acid: LD50 (mammal, skin)=1060mg/kg; LD50 (rabbit, skin)=1060 mg/kg; LC50(inhalation, mouse)=5620 ppm/1H; LC50(inhalation, mouse)=5620 mg/m<sup>3</sup>/1H  
Ethylene Glycol: LD50(oral, cat)=1650 mg/kg, LD50 (oral, cat)=2000 mg/kg

**Carcinogenicity (NTP, IARC, OSHA):** No reagents listed as carcinogens

### 12. Ecological Information

**Ecotoxicity:** Acute Aquatic Effects Data for 100% Glacial Acetic Acid 96 h LC-50 (fathead minnow): > 100mg/L 48 h LC-50 (golden orfe): 410 mg/L 48 h LC-50 (mosquito fish): 251 mg/L 96 h LC-50 (daphnid): > 100 mg/L

**Persistence and degradability:** No long term effects known.

**Bioaccumulation Potential (octanol-water partition coefficient, BCF):** Oxygen Demand Data for 100% Glacial Acetic Acid BOD-5: 340-880 mg/g BOD-20: 900 mg/g COD: 1,030 mg/g

## Harris Hematoxylin

**Mobility in the soil:** N/A

**Adverse Environmental Effects:** N/A

### 13. Disposal Considerations

**Recommended Disposal Containers:** Check with your local waste authorities\*

**Recommended Disposal Methods:** Do not dispose of in drains, check with your local waste authorities.\*

**Physical/Chemical Properties affecting Disposal:** See section 2 and section 9 applicable information.\*

**Special Precautions for Landfill and Incineration Activities:** Check with your local waste authorities.\*

**Waste Stream:** Consult your local or regional authorities.\*

### 14. Transport Information

<b>14.1 UN Number DOT, IATA,IMDG, ADR</b>	<b>DOT:</b> Not Regulated for Transport by DOT. <b>ADR/RID:</b> Not Regulated for Transport by ADR/RID. <b>IMDG:</b> Not Regulated for Transport by IMDG. <b>IATA:</b> Not Regulated for Transport by IATA.
<b>14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR</b>	<b>ADR/RID:</b> Not Regulated for Transport by ADR/RID. <b>IMDG:</b> Not Regulated for Transport by IMDG. <b>IATA:</b> Not Regulated for Transport by IATA. <b>DOT:</b> Not Regulated for Transport by DOT.
<b>14.3 Transport Hazard Class(es)</b>	<b>DOT:</b> Not Regulated for Transport by DOT. <b>IATA:</b> Not Regulated for Transport by IATA. <b>IMDG:</b> Not Regulated for Transport by IMDG. <b>ADR/RID:</b> Not Regulated for Transport by ADR/RID.
<b>14.4 Packing Group DOT, IATA,IMDG, ADR</b>	<b>ADR/RID:</b> Not Regulated for Transport by ADR/RID. <b>IMDG:</b> Not Regulated for Transport by IMDG. <b>IATA:</b> Not Regulated for Transport by IATA. <b>DOT:</b> Not Regulated for Transport by DOT.
<b>14.5 Environmental Hazards</b>	<b>Marine Pollutant:</b> No
<b>14.6 Special Precautions for User</b>	Not applicable.

### 15. Regulatory Information

**OSHA:**

**DOT:**

**EPA:**

**CPSC:**

## Harris Hematoxylin


### 16. Other Information

Revision Date: 2019-01-11

#### NFPA

Health	1
Fire Hazard	1
Reactivity	1
Specific Hazard	

National Fire Protection Association (USA) NFPA



#### HMIS

Health	1
Flammability	1
Physical Hazard	1
Personal Protection	

Hazardous Material Information System HMIS

Health	1
Flammability	1
Physical Hazard	1
Personal Protection	

#### Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.