

## 1. Identification

### Product identifier

**Product Name** PureAffin® i

### Other means of identification

**Catalogue Number** EEPAR1

**Synonyms** None

### Recommended use of the chemical and restrictions on use

**Recommended use** Laboratory chemicals

**Restrictions on use** No information available

### Details of the supplier of the safety data sheet

#### **Company**

Cancer Diagnostics, Inc.  
4300 Emperor Blvd. #400  
Durham, NC 27703  
1-877-846-5393

### Emergency telephone number

**Emergency Telephone** Chemtrec US: (800) 424-9300  
Infotrac International  
1-352-323-3500

## 2. Hazard(s) identification

### Classification

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

### Hazards not otherwise classified (HNOC)

Not applicable

### Label elements

#### Hazard statements

Not a hazardous substance or mixture according to the Globally Harmonized System (GHS)

The product contains no substances which at their given concentration, are considered to be hazardous to health.

**Appearance** white**Physical state** Solid**Odor** Slight**Other information**

May be harmful in contact with skin

**Unknown acute toxicity** 99.12 % of the mixture consists of ingredient(s) of unknown toxicity  
 0.01 % of the mixture consists of ingredient(s) of unknown acute oral toxicity  
 0.01 % of the mixture consists of ingredient(s) of unknown acute dermal toxicity  
 99.12 % of the mixture consists of ingredient(s) of unknown acute inhalation toxicity (gas)  
 99.12 % of the mixture consists of ingredient(s) of unknown acute inhalation toxicity (vapor)  
 99.12 % of the mixture consists of ingredient(s) of unknown acute inhalation toxicity (dust/mist)

**3. Composition/information on ingredients****Substance**

Not applicable.

**Mixture**

Chemical name	CAS No.	Weight-%
Paraffin waxes and Hydrocarbon waxes	8002-74-2	>99
Polyisobutylene	9003-27-4	<1
Dimethyl sulfoxide	67-68-5	<1
2,6-Di-tert-butyl-p-cresol	128-37-0	<.5

**4. First-aid measures****Description of first aid measures**

**Inhalation** Remove to fresh air.

**Eye contact** Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a physician.

**Skin contact** Wash skin with soap and water.

**Ingestion** Clean mouth with water and drink afterwards plenty of water.

**Most important symptoms and effects, both acute and delayed**

**Symptoms** No information available.

**Indication of any immediate medical attention and special treatment needed**

**Note to physicians** Treat symptomatically.

**5. Fire-fighting measures**

**Suitable Extinguishing Media** Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

**Unsuitable extinguishing media** CAUTION: Use of water spray when fighting fire may be inefficient.

**Specific hazards arising from the chemical** No information available.

**Explosion data**

**Sensitivity to mechanical impact** None.  
**Sensitivity to static discharge** None.

**Special protective equipment for fire-fighters** Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

**6. Accidental release measures**

**Personal precautions, protective equipment and emergency procedures**

**Personal precautions** Ensure adequate ventilation.

**Methods and material for containment and cleaning up**

**Methods for containment** Prevent further leakage or spillage if safe to do so.

**Methods for cleaning up** Pick up and transfer to properly labeled containers.

**7. Handling and storage**

**Precautions for safe handling**

**Advice on safe handling** Handle in accordance with good industrial hygiene and safety practice.

**Conditions for safe storage, including any incompatibilities**

**Storage Conditions** Keep containers tightly closed in a dry, cool and well-ventilated place.

**8. Exposure controls/personal protection**

**Control parameters**

**Exposure Limits**

Chemical name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Paraffin waxes and Hydrocarbon waxes 8002-74-2	TWA: 2 mg/m <sup>3</sup> fume	(vacated) TWA: 2 mg/m <sup>3</sup>	TWA: 2 mg/m <sup>3</sup> fume
2,6-Di-tert-butyl-p-cresol 128-37-0	TWA: 2 mg/m <sup>3</sup> inhalable fraction and vapor	(vacated) TWA: 10 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup>

**Appropriate engineering controls**

**Engineering controls** Showers  
 Eyewash stations  
 Ventilation systems.

**Individual protection measures, such as personal protective equipment**

**Eye/face protection** No special protective equipment required.

**Skin and body protection** No special protective equipment required.

**Respiratory protection** No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.

**General hygiene considerations** Handle in accordance with good industrial hygiene and safety practice.

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state	Solid
Appearance	white
Color	No information available
Odor	Slight
Odor threshold	No information available

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
pH	No data available	None known
Melting point / freezing point	55-57 °C / °F	
Boiling point / boiling range	315.5 °C / 599.9 °F	
Flash point	204.39 °C / 399.9 °F	
Evaporation rate	No data available	None known
Flammability (solid, gas)	No data available	None known
Flammability Limit in Air		None known
Upper flammability or explosive limits	No data available	
Lower flammability or explosive limits	No data available	
Vapor pressure	No data available	None known
Vapor density	Not applicable	
Relative density	No data available	None known
Water solubility	No data available	None known
Solubility in other solvents	No data available	None known
Partition coefficient	No data available	None known
Autoignition temperature	No data available	None known
Decomposition temperature	No data available	None known
Kinematic viscosity	No data available	None known
Dynamic viscosity	No data available	

### Other information

Explosive properties	No information available
Oxidizing properties	No information available
Softening point	No information available
Molecular weight	No information available
VOC Content (%)	0.31
Liquid Density	No information available
Bulk density	No information available

## 10. Stability and reactivity

Reactivity	No information available.
Chemical stability	Stable under normal conditions.
Possibility of hazardous reactions	None under normal processing.
Conditions to avoid	None known based on information supplied.
Incompatible materials	None known based on information supplied.
Hazardous decomposition products	None known based on information supplied.

## 11. Toxicological information

### Information on likely routes of exposure

**Product Information**

**Inhalation** Specific test data for the substance or mixture is not available.  
**Eye contact** Specific test data for the substance or mixture is not available.  
**Skin contact** May be harmful in contact with skin.  
**Ingestion** Specific test data for the substance or mixture is not available.

**Symptoms related to the physical, chemical and toxicological characteristics**

**Symptoms** No information available.

**Acute toxicity**

**Numerical measures of toxicity**

The following values are calculated based on chapter 3.1 of the GHS document .

**ATEmix (oral)** 5,049.94 mg/kg  
**ATEmix (dermal)** 3,636.00 mg/kg

**Unknown acute toxicity** 99.12 % of the mixture consists of ingredient(s) of unknown toxicity  
 0.01 % of the mixture consists of ingredient(s) of unknown acute oral toxicity  
 0.01 % of the mixture consists of ingredient(s) of unknown acute dermal toxicity  
 99.12 % of the mixture consists of ingredient(s) of unknown acute inhalation toxicity (gas)  
 99.12 % of the mixture consists of ingredient(s) of unknown acute inhalation toxicity (vapor)  
 99.12 % of the mixture consists of ingredient(s) of unknown acute inhalation toxicity (dust/mist)

**Component Information**

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Paraffin waxes and Hydrocarbon waxes 8002-74-2	> 5000 mg/kg ( Rat )	> 3600 mg/kg ( Rabbit )	-
Dimethyl sulfoxide 67-68-5	= 28300 mg/kg ( Rat ) = 14500 mg/kg ( Rat )	= 40 g/kg ( Rat )	> 5.33 mg/L ( Rat ) 4 h
2,6-Di-tert-butyl-p-cresol 128-37-0	> 2930 mg/kg ( Rat )	> 2000 mg/kg ( Rat )	-

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

**Skin corrosion/irritation** No information available.  
**Serious eye damage/eye irritation** No information available.  
**Respiratory or skin sensitization** No information available.  
**Germ cell mutagenicity** No information available.  
**Carcinogenicity** No information available.

Chemical name	ACGIH	IARC	NTP	OSHA
2,6-Di-tert-butyl-p-cresol 128-37-0	-	Group 3	-	-

**Reproductive toxicity** No information available.  
**STOT - single exposure** No information available.

<b>STOT - repeated exposure</b>	No information available.
<b>Target organ effects</b>	Respiratory system, Eyes, Skin.
<b>Aspiration hazard</b>	No information available.
<b>Other adverse effects</b>	No information available.
<b>Interactive effects</b>	No information available.

## 12. Ecological information

### Ecotoxicity

Chemical name	Algae/aquatic plants	Fish	Toxicity to microorganisms	Crustacea
Dimethyl sulfoxide 67-68-5	EC50: 12350 - 25500mg/L (96h, Skeletonema costatum)	LC50: >40g/L (96h, Lepomis macrochirus) LC50: =34000mg/L (96h, Pimephales promelas) LC50: 33 - 37g/L (96h, Oncorhynchus mykiss) LC50: =41.7g/L (96h, Cyprinus carpio)	-	EC50: =7000mg/L (24h, Daphnia species)
2,6-Di-tert-butyl-p-cresol 128-37-0	EC50: >0.42mg/L (72h, Desmodemus subspicatus) EC50: =6mg/L (72h, Pseudokirchneriella subcapitata)	LC50: =5mg/L (48h, Oryzias latipes)	-	-

**Persistence and degradability** No information available.

**Bioaccumulation** There is no data for this product.

Chemical name	Partition coefficient
Dimethyl sulfoxide 67-68-5	-2.03
2,6-Di-tert-butyl-p-cresol 128-37-0	4.17

### Mobility

**Other adverse effects** No information available.

## 13. Disposal considerations

### Waste treatment methods

**Waste from residues/unused products** Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.

**Contaminated packaging** Do not reuse empty containers.

## 14. Transport information

<b>DOT</b>	Not regulated
<b>TDG</b>	Not regulated
<b>IATA</b>	Not regulated
<b>IMDG</b>	Not regulated

## 15. Regulatory information

### International Inventories

<b>TSCA</b>	Contact supplier for inventory compliance status.
<b>DSL/NDSL</b>	Contact supplier for inventory compliance status.
<b>EINECS/ELINCS</b>	Contact supplier for inventory compliance status.
<b>ENCS</b>	Contact supplier for inventory compliance status.
<b>IECSC</b>	Contact supplier for inventory compliance status.
<b>KECL</b>	Contact supplier for inventory compliance status.
<b>PICCS</b>	Contact supplier for inventory compliance status.
<b>AICS</b>	Contact supplier for inventory compliance status.

### **Legend:**

- TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory
- DSL/NDSL** - Canadian Domestic Substances List/Non-Domestic Substances List
- EINECS/ELINCS** - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances
- ENCS** - Japan Existing and New Chemical Substances
- IECSC** - China Inventory of Existing Chemical Substances
- KECL** - Korean Existing and Evaluated Chemical Substances
- PICCS** - Philippines Inventory of Chemicals and Chemical Substances
- AICS** - Australian Inventory of Chemical Substances

### US Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

#### **SARA 311/312 Hazard Categories**

Should this product meet EPCRA 311/312 Tier reporting criteria at 40 CFR 370, refer to Section 2 of this SDS for appropriate classifications. Under the amended regulations at 40 CFR 370, EPCRA 311/312 Tier II reporting for the 2017 calendar year will need to be consistent with updated hazard classifications.

#### **CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

#### **CERCLA**

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material.

### US State Regulations

#### **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

**U.S. State Right-to-Know Regulations**

**US State Regulations**

Chemical name	New Jersey	Massachusetts	Pennsylvania
Paraffin waxes and Hydrocarbon waxes 8002-74-2	X	-	X
Dimethyl sulfoxide 67-68-5	X	-	-
2,6-Di-tert-butyl-p-cresol 128-37-0	X	-	X

**U.S. EPA Label Information**

**EPA Pesticide Registration Number** Not applicable

**16. Other information**

<b>NFPA</b>	<b>Health hazards</b> 0	<b>Flammability</b> 1	<b>Instability</b> 0	<b>Physical and chemical properties</b> -
<b>HMIS</b>	<b>Health hazards</b> 0	<b>Flammability</b> 1	<b>Physical hazards</b> 0	<b>Personal protection</b> X

**Key or legend to abbreviations and acronyms used in the safety data sheet**

**Legend Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION**

TWA	TWA (time-weighted average)	STEL	STEL (Short Term Exposure Limit)
Ceiling	Maximum limit value	*	Skin designation

**Key literature references and sources for data used to compile the SDS**

Agency for Toxic Substances and Disease Registry (ATSDR)  
 U.S. Environmental Protection Agency ChemView Database  
 European Food Safety Authority (EFSA)  
 EPA (Environmental Protection Agency)  
 Acute Exposure Guideline Level(s) (AEGl(s))  
 U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act  
 U.S. Environmental Protection Agency High Production Volume Chemicals  
 Food Research Journal  
 Hazardous Substance Database  
 International Uniform Chemical Information Database (IUCLID)  
 Japan GHS Classification  
 Australia National Industrial Chemicals Notification and Assessment Scheme (NICNAS)  
 NIOSH (National Institute for Occupational Safety and Health)  
 National Library of Medicine's ChemID Plus (NLM CIP)  
 National Library of Medicine's PubMed database (NLM PUBMED)  
 National Toxicology Program (NTP)  
 New Zealand's Chemical Classification and Information Database (CCID)  
 Organization for Economic Co-operation and Development Environment, Health, and Safety Publications  
 Organization for Economic Co-operation and Development High Production Volume Chemicals Program  
 Organization for Economic Co-operation and Development Screening Information Data Set  
 RTECS (Registry of Toxic Effects of Chemical Substances)  
 World Health Organization

<b>Issuing Date</b>	20-Feb-2014
<b>Revision date</b>	08-Oct-2020
<b>Revision Note</b>	SDS sections updated.



**Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

**End of Safety Data Sheet**

Catalogue Number

EEPAR1