



1 Well labelled and easy to read identification of biopsy origin

2 Each biopsy is isolated from each other by a divider

3 Molded ridges around lid make it easy to open and close

4 Fast thread design allows lid to be removed with only 1 ½ turn

5 Specially designed O-ring for a more secure and positive leakproof seal

6 Warning label has space for patient identification

LEFT SIDE

RIGHT SIDE

PROSTATE BIOPSY

Patient: _____

ID #: _____ Date: _____

Dr.: _____ Time: _____

IVD

Patent Pending www.simport.com Made in Canada

CE