




Section 1. Identification

1.1 Product Identifier	Product Name: Sodium Metabisulfite Solution (5%) Product Number: SSC1129	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 Classification: Regulation (EC) No 1272/2008 (CLP)	Serious eye damage (Category 1) – H318	
Classification: OSHA 29 CFR 1910.1200 HCS	Serious eye damage (Category 1) – H318	
Classification: Directive 1999/45/EC		Xn R22, R31, R41, S2, S26, S39, S46 EN: Harmful DA: Sundhedsskadelig DE: Gesundheitsschädlich EL: Επιβλαβές ES: Nocivo FR: Nocif IT: Nocivo NL: Schadelijk PT: Nocivo SV: Hälsoskadlig For full text of R- & S- phrases: see Section 15.


2.2 Label Elements	Pictogram(s):  Signal word: Danger Hazard statement(s): H318 – Causes serious eye damage Precautionary statement(s): P302+ P352 – IF ON SKIN: Wash with plenty of soap and water P305 + P351 – IF IN EYES: Rinse cautiously with water for several minutes. P337 + P313 – If eye irritation persists: Get medical advice/attention								
NFPA Scale: 0-4									
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr> <td>HEALTH</td> <td>2</td> </tr> <tr> <td>FLAMMABILITY</td> <td>0</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>C</td> </tr> </table>	HEALTH	2	FLAMMABILITY	0	PHYSICAL HAZARD	0	PERSONAL PROTECTION	C
HEALTH	2								
FLAMMABILITY	0								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	C								
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.								

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Sodium Metabisulfite	7681-57-4	231-673-0	 Danger. 1 H318 Causes serious eye damage. Warning. 4 H302 Harmful if swallowed.	≤ 5

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water. Get immediate medical attention if irritation continues.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get medical attention if symptoms develop.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	Releases toxic sulfur oxides under fire conditions.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Sodium Metabisulfite: NIOSH REL: TWA 5 mg/m ³ ACGIH TLV: TWA 5 mg/m ³
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor.

	Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Transparent.
Odor	Strong sulfur odor.
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Oxidizing agents.
10.6 Hazardous Decomposition Materials	Toxic sulfur oxides under fire conditions.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. May be irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Corrosive to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged unless thoroughly diluted. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	None.

Section 14. Transport Information

14.1 UN Number DOT, IATA, IMDG, ADR	DOT: Not dangerous goods. ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods.
14.2 UN Proper Shipping Name DOT, IATA, IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods. IATA: Not dangerous goods. IMDG: Not dangerous goods.

	ADR/RID: Not dangerous goods.
14.4 Packing Group DOT, IATA, IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	<p>R22 EN: Harmful if swallowed. DA: Farlig ved indtagelse. DE: Gesundheitsschädlich beim Verschlucken. EL: Επιβλαβές σε περίπτωση καταπόσεως. ES: Nocivo por ingestión. FR: Nocif en cas d'ingestion. IT: Nocivo per ingestione. NL: Schadelijk bij opname door de mond. PT: Nocivo por ingestão. SV: Farligt vid förtäring.</p> <p>R31 EN: Contact with acids liberates toxic gas. DA: Udvikler giftig gas ved kontakt med syre. DE: Entwickelt bei Berührung mit Säure giftige Gase. EL: Σε επαφή με οξέα ελευθερώνονται τοξικά αέρια. ES: Σε επαφή με οξέα ελευθερώνονται τοξικά αέρια. FR: Au contact d'un acide, dégage un gaz toxique. IT: A contatto con acidi libera gas tossico. NL: Vormt vergiftige gassen in contact met zuren. PT: Em contacto com ácidos liberta gases tóxicos. SV: Utvecklar giftig gas vid kontakt med syra.</p> <p>R41 EN: Risk of serious damage to eyes. DA: Risiko for alvorlig øjenskade. DE: Gefahr ernster Augenschäden. EL: Κίνδυνος σοβαρών οφθαλμικών βλαβών. ES: Riesgo de lesiones oculares graves. FR: Risque de lésions oculaires graves. IT: Rischio di gravi lesioni oculari. NL: Gevaar voor ernstig oogletsel. PT: Risco de lesões oculares graves. SV: Risk för allvarliga ögonskador.</p>
S phrases	<p>S2 EN: Keep out of the reach of children.</p>

	<p>DA: Opbevares utilgængeligt for børn. DE: Darf nicht in die Hände von Kindern gelangen. EL: Μακριά από παιδιά. ES: Manténgase fuera del alcance de los niños. FR: Conserver hors de la portée des enfants. IT: Conservare fuori della portata dei bambini. NL: Buiten bereik van kinderen bewaren. PT: Manter fora do alcance das crianças. SV: Förvaras oåtkomligt för barn.</p> <p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετέ τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen. PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista. SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S39 EN: Wear eye/face protection. DA: Brug beskyttelsesbriller/ansigtsskærm under arbejdet. DE: Schutzbrille/Gesichtsschutz tragen. EL: Χρησιμοποιείτε συσκευή προστασίας ματιών/προσώπου. ES: Úsese protección para los ojos/la cara. FR: Porter un appareil de protection des yeux/du visage. IT: Proteggersi gli occhi/la faccia. NL: Een bescherming voor de ogen/voor het gezicht dragen. PT: Usar um equipamento protector para os olhos /face. SV: Använd skyddsglasögon eller ansiktsskydd.</p> <p>S46 EN: If swallowed, seek medical advice immediately and show this container or label. DA: Ved indtagelse, kontakt odümgående læge og vis denne beholder eller etiket. DE: Bei Verschlucken sofort ärztlichen Rat einholen und Verpackung odor Etikett vorzeigen. EL: Σε περίπτωση κατάποσης ζητήστε αμέσως ιατρική συμβουλή και δείξτε αυτό το δοχείο ή την ετικέτα. ES: En caso de ingestión, acúdase inmediatamente al medico y muéstrele la etiqueta o el envase. FR: En cas d'ingestion, consulter immédiatement un médecin et lui montrer l'emballage ou l'étiquette. IT: In caso d'ingestione consultare immediatamente il medico e mostrargli il contenitore o l'etichetta. NL: In geval van inslikken onmiddellijk een arts raadplegen en verpakking of etiket tonen. PT: Em caso de ingestão, consultar imediatamente o médico e mostrar-lhe a embalagem ou o rótulo. SV: Vid förtäring kontakta genast läkare och visa denna förpackning eller etiketten.</p>
--	---

Section 16. Other Information

Revision Date: 2019-01-11



Safety Data Sheet

Emergency 1-800-424-9300
Contact (intl.) 1-352-323-3500

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.