SAFETY DATA SHEET



Revision Date 2022-06-08 Version 2

1. IDENTIFICATION

Product Name Ribbon Tissue Section Adhesive

Product Code P1000RD, PCM1005

Recommended Use For laboratory use only.

Company Cancer Diagnostics, Inc.

116 Page Point Circle Durham, NC 27703 Tel. (877) 846-5393

Emergency Telephone Call CHEMTREC 1-800-424-9300

2. HAZARDS IDENTIFICATION

Classification

OSHA Regulatory Status

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Label elements

Although this product does not normally require specific hazard precautions, chemical users should always take care to minimize personnel exposure and workplace contamination. Store tightly closed at 15-30C.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS No.	Weight-%
Water	7732-18-5	>95
Gelatin	9000-70-8	<2
Chromium potassium sulfate dodecahydrate	7788-99-0	<2

4. FIRST AID MEASURES

Description of first aid measures

General advice This chemical is not expected to produce any significant adverse health effects.

Eye contact Flush eyes with plenty of water, removing contact lenses if present. Get medical attention if

irritation develops.

Skin contact Wash thoroughly with soap and water while removing contaminated garments. Get medical

attention if irritation develops.

Inhalation Remove to fresh air. Get medical attention for any breathing difficulty.

Ingestion Rinse mouth and drink several glasses of water. Contact a physician or poison control

center if symptoms develop.

5. FIRE-FIGHTING MEASURES

Suitable extinguishing media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment

Specific hazards arising from the chemical

No information available.

Hazardous combustion products Carbon oxides. Sulphur oxides. Oxides of potassium.

Protective equipment and precautions for firefighters

Firefighters should wear self-contained breathing apparatus with full facepiece operated in pressure-demand or other positive pressure mode.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal precautions Use personal protective equipment as required. Ensure adequate ventilation, especially in

confined areas.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Absorb spill with inert material, scoop up and containerize for disposal.

7. HANDLING AND STORAGE

Handle in accordance with good industrial hygiene and safety practice.

Storage Conditions Keep containers tightly closed in a dry, cool and well-ventilated place.

Incompatible materials None known based on information supplied.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Occupational exposure limits
Appropriate engineering controls

Not applicable

Engineering Controls Emergency showers, eyewash stations, ventilation systems.

General Hygiene Considerations Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical state Liquid

AppearanceClear to slightly grey-greenOdorNo information availableOdor thresholdNo information available

pH No information available
Melting point / freezing point
Boiling point / boiling range
Flash point No information available
No information available
Evaporation rate No information available

Ribbon Tissue Section Adhesive

Flammability (solid, gas) No information available

Flammability Limit in Air

Upper flammability limit: No information available No information available Lower flammability limit: No information available Vapor pressure No information available Vapor density Relative density No information available Miscible with water Water solubility No information available Solubility in other solvents No information available **Partition coefficient Autoignition temperature** No information available **Decomposition temperature** No information available Kinematic viscosity No information available

10. STABILITY AND REACTIVITY

Stability Stable under recommended storage conditions.

Possibility of Hazardous Reactions No information available.

Conditions to avoid Extremes of temperature and direct sunlight

Incompatible materials None known based on information supplied.

Hazardous Decomposition Products Carbon oxides. Sulphur oxides. Oxides of potassium.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Symptoms No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation
Serious eye damage/eye irritation
Irritation
Sensitization
Sensitization
Germ cell mutagenicity

No information available.
No information available.
No information available.
No information available.

Carcinogenicity The table below indicates whether each agency has listed any ingredient as a carcinogen

IARC (International Agency for Research on Cancer) Group 3 - Not classifiable as to carcinogenicity in humans

12. ECOLOGICAL INFORMATION

Ecotoxicity

None known

Persistence and degradability

No information available.

Bioaccumulation

No information available.

13. DISPOSAL CONSIDERATIONS

Waste treatment methods

Disposal of wastes

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Chemical Name	California Hazardous Waste Status	
Chromium potassium sulfate dodecahydrate	Toxic	
7788-99-0	Corrosive	
	Ignitable	

14. TRANSPORT INFORMATION

Note: Transportation information is provided as a general reference only and may not be applicable in all situations. This information applies to non-bulk shipments only. Per 49 CFR §173.22, it is the shipper's responsibility to ensure that all materials are properly packaged, classified and labeled prior to shipment.

DOT Not regulated

IATA Not regulated

15. REGULATORY INFORMATION

No information available.

16. OTHER INFORMATION

Prepared By CDI Regulatory Affairs (Email: compliance@cancerdiagnostics.com)

Revision Date 2022-06-08

Disclaimer

Cancer Diagnostics, Inc. makes no representation as to the comprehensiveness or accuracy of this document. Individuals using this information, or the product to which it refers, must exercise their independent judgment in determining all appropriateness for a particular purpose. Accordingly, Cancer Diagnostics, Inc. will not be responsible for damages of any kind or nature resulting from the use of this information or the corresponding product. No representations or warranties of any kind or nature, including but not limited to: express warranties, implied warranties or merchantability, or warranties of fitness for a particular purpose, are made hereunder with respect to the information set forth herein or to the product to which the information refers.

End of Safety Data Sheet