

Clear-Away Citrus™

1. Identification

Product Name: Clear-Away Citrus™

Item #: FX2179, FX2180, FX2181

Synonyms: D-limonene; Hemo-De; Carvene; (+) -4- Isopropenyl -1-methylcyclohexene; (+) -R-Limonene; (R) - 1- Methyl -4- (1-methylethenyl) cyclohexene; D- (+) - Limonene; d-p-Mentha-1, 8-diene; optically active terpene; (The racemic mixture of d-limonene and l-limonene isomers is known as dipentene).

Recommended Use: Clearing Agent

Manufacturer:

 Cancer Diagnostics, Inc.
 116 Page Point Circle
 Durham, NC 27703
 1-877-846-5393

Restrictions on Use: N/A

In Case of Emergency:

 Chemtrec US 1-800-424-9300
 Chemtrec International 703-527-3887

2. Hazards Identification

OSHA Hazard Classification(s):

 Skin Irritation - Category 2
 Eye Irritation - Category 2A
 Sensitization - Skin - Category 1B
 Flammable Liquids - Category 3

Signal Word: Warning

Hazard Statement(s): Causes skin irritation. Causes serious eye irritation. May cause an allergic skin reaction. Flammable liquid and vapor.

Pictogram(s):

Precautionary Statement(s): Prevention: Wash body thoroughly after handling. Wear protective gloves. Wear eye protection, face protection. Avoid breathing dust, vapors. Contaminated work clothing must not be allowed out of the workplace. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

Response: If on skin: Wash with plenty of water. Specific treatment (see first aid section on this label). If skin irritation or rash occurs: Get medical attention. Take off all contaminated clothing and wash it before reuse. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing If eye irritation persists: Get medical attention. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. In case of fire: Use water, dry chemical, CO2 or foam to extinguish.

Storage: Store in a well-ventilated place. Keep cool.

Disposal: Dispose of contents/container in accordance with local regulations.

Descriptions of Hazards not otherwise classified: N/A

Percent of mixture with unknown acute toxicity: N/A

3. Composition and Information on Ingredients

Chemical Name	Common Name	CAS #	Concentration %
D-limonene		5989-27-5	Greater than 95%

4. First Aid Measures

Eye Contact: If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Clear-Away Citrus™

Continue rinsing. If eye irritation persists: Get medical advice/attention.

Skin Contact: If on skin (or hair): Take off immediately all contaminated clothing and wash before reuse. Wash with plenty of water. If skin irritation occurs: Get medical advice/attention.

Inhalation: Remove to fresh air; give artificial respiration if breathing has stopped. Get medical advice/attention if you feel unwell.

Ingestion: Do not induce vomiting. If vomiting occurs spontaneously, keep head below hips to prevent aspiration of liquid into the lungs. Get medical attention.

Symptoms: Irritation eyes, nose, throat; headache, dizziness. May cause skin sensitization.

Recommendations for immediate medical care/special treatment: Get medical advice/attention if you feel unwell.

5. Fire- Fighting Measures

Extinguishing Media: Dry chemical, carbon dioxide, alcohol foam, water. Use water spray to cool fire-exposed containers and disperse vapors.

Fire Hazards (Chemical): OSHA classified Flammable Liquid

Special Protective Equipment: Fire fighters should use self-contained breathing apparatus and protective clothing.

Precautions for Firefighters: Carbon monoxide and unidentified organic compounds may be formed during combustion. Vapors can travel distances to ignition source and flash back. Cool fire exposed containers with water. Fine mist or spray may be flammable at temperatures below the flash point. When heated above the flash point this material emits flammable vapors which, when mixed with air, can burn or be explosive.

6. Accidental Release Measures

Emergency Procedures: Evacuate the area of all unnecessary personnel. Wear suitable protective equipment. Eliminate all sources of ignition and provide ventilation.

Protective Equipment: See section 8

Environmental Precautions: Prevent release to the environment by using barriers.

Containment and Clean-Up Procedures: Use barriers to prevent spreading. Collect spill in container. Call waste authorities.

7. Handling and Storage

Handling: Do not breathe vapors. Do not eat, drink or smoke when using this product. Keep away from heat, sparks, open flames, hot surfaces. Keep container tightly closed. Ground/Bond container and receiving equipment. Use explosion-proof electrical, ventilating, and lighting equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

Storage: Store in a well-ventilated place. Keep cool. Store away from sources of ignition, open flame and heat. Keep container tightly sealed when not in use. Separate from oxidizing materials. Partially filled containers should be blanketed with nitrogen.

8. Exposure Controls/Personal Protection

OSHA Permissible Exposure Limits (PELs):

Reagent	CAS #	OSHA PEL TWA
N/A		

ACGIH Threshold Limit Values (TLVs):

Reagent	CAS #	ACGIH PEL TLV	ACGIH STEL
N/A			

Engineering Controls: Use in a well ventilated area to prevent exposure. Maintain eyewash fountain and quick-drench facilities in work areas.

Personal Protective Measures: Wear gloves, lab coat, eye protection and impervious footwear. Contact lenses should not be worn when working with this material.

Special PPE Requirements: If ventilation hood not available wear respirator.

9. Physical and Chemical Properties Section

Appearance: Almost Colorless, Liquid

Clear-Away Citrus™

Molecular Weight: 136.2
Molecular Formula: C₁₀H₁₆
pH: N/A
Boiling Point and Boiling Range: 177°C
Melting Point/Freezing Point: -96.9°C
Flash Point: 43°C
Specific Gravity/Relative Density: 0.8400
Odor: Citrus-like aroma
Odor Threshold: N/A
Color: Almost Colorless
Flammability (solid/gas): N/A
Vapor Density: N/A
Upper/Lower flammability or explosive limits: LEL: 0.70 vol%, UEL: 6.10 vol%
Vapor Pressure: N/A
Evaporation Rate: N/A
Partition Coefficient: n-octanol/water: N/A
Viscosity: N/A
Auto-ignition temperature: 255°C
Solubility: Insoluble in Water.
Decomposition Temperature: N/A

10. Stability and Reactivity

Reactivity:
Chemical Stability: Stable
Conditions of Stability/Instability: Stable under normal temperatures and pressures. Oxidizes to a film in air, oxidation behavior similar to that of rubber or drying oils. Avoid heat, flame, and other sources of ignition and prolonged air exposure.
Stabilizers needed: None
Safety issue indicated by appearance change: N/A
Other: N/A
Hazardous Reactions: N/A
Hazardous Polymerization: Does not occur
Conditions to avoid: N/A
Classes of Incompatible Materials: Oxidizers, Strong Acids, Strong Bases
Hazardous Decomposition Products: Thermal-oxidation degradation can produce oxides of carbon. Toxic gases and vapors (i.e. Carbon monoxide) may be released in a fire.

11. Toxicological Information**Likely Routes of Exposure**

Eyes: Irritation.

Skin: Irritation. May cause skin sensitization, an allergic reaction, which becomes evident upon re-exposure to this material. Prolonged and/ or repeated contact may cause defatting of the skin and dermatitis.

Inhalation: Dizziness, headache. Irritation of respiratory tract.

Ingestion: Nausea. May cause digestive tract disturbances.

Signs or Symptoms of Exposure: Nausea. Skin Irritation.

Effects from short term exposure (delayed, immediate, chronic): Irritation to the eyes, nose, throat; headache, dizziness, nausea.

Acute Toxicity (Numerical Measures): N/A

Carcinogenicity (NTP, IARC, OSHA): IARC Group 3 "Not classifiable as to its carcinogenicity to humans"

Clear-Away Citrus™**12. Ecological Information****Ecotoxicity:** N/A**Persistence and degradability:** N/A**Bioaccumulation Potential (octanol-water partition coefficient, BCF):** N/A**Mobility in the soil:** N/A**Adverse Environmental Effects:** N/A**13. Disposal Considerations****Recommended Disposal Containers:** Check with your local waste authorities***Recommended Disposal Methods:** Do not dispose of in drains, check with your local waste authorities.***Physical/Chemical Properties affecting Disposal:** See section 2 and section 9 applicable information.***Special Precautions for Landfill and Incineration Activities:** Check with your local waste authorities.***Waste Stream:** Consult your local or regional authorities.***14. Transport Information****UN Number:** UN2319**UN Proper Shipping Name:** Terpene Hydrocarbons, n.o.s., (D-Limonene)**Transport Hazard Class(es):** 3**Packing Group Number:** III**Environmental Hazards (IMDG code):****Marine Pollutant:** No**Transport in Bulk (IBC Code):** N/A**Special Transport Precautions:** N/A**15. Regulatory Information****OSHA:****DOT:****EPA:****CPSC:**

Clear-Away Citrus™

16. Other Information


Revision Date: 2019-01-11

NFPA

Health	1
Fire Hazard	3
Reactivity	0
Specific Hazard	

National Fire Protection Association (USA) NFPA

Fire Hazard



HMIS

Health	1
Flammability	3
Physical Hazard	0
Personal Protection	

Hazardous Material Information System HMIS

Health	1
Flammability	3
Physical Hazard	0
Personal Protection	

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.