



# Pap Spray™

ALCOHOLIC PAP SMEAR FIXATIVE

#FX1086-2 2OZ SPRAY PUMP x 12/CS

#FX1086 5OZ SPRAY PUMP x 12/CS

Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_

Pap Spray™ is a polyethylene glycol formulation for rapid, high quality fixation of cytology cell spreads. A convenient fine-mist sprayer covers cells with a layer that protects cell morphology yet is water and alcohol soluble.



## DANGER

### HAZARD STATEMENTS

Causes skin and serious eye irritation. Causes damage to organs (respiratory system). May cause damage to organs (respiratory system, central nervous system, liver, blood) through prolonged or repeated exposure. Highly flammable liquid and vapor.

### PRECAUTIONARY STATEMENTS • PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Do not breathe dust, vapors. Do not eat, drink, or smoke when using this product. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

### PRECAUTIONARY STATEMENTS • RESPONSE

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. If eye irritation persists: get medical attention. If exposed or concerned: call a doctor. Call a doctor if you feel unwell.

**IF ON SKIN:** Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

**IN CASE OF FIRE:** Use CO<sub>2</sub>, dry chemical, or foam for extinction.

### STORAGE & DISPOSAL

Store locked up. Store in a well-ventilated place. Keep cool. Dispose of contents/container to an approved waste disposal plant.

**NOTICE:** Read Safety Data Sheet for additional information.

**EMERGENCY CONTACT:** Chemtec 1-800-424-9300

For In Vitro Diagnostics Use



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)

Made in USA