



Zorbentz™ Lint-Free
Disposable Cloths
Datasheet

“The Cure for
the Common Lab”

Properties					
Item Code	Color	Ultimate Tensile Strength	Thickness - Caliper (single folded sheet) 0.01mm	Thickness - Caliper (50 sheets, banded) 0.01mm	Total Absorbant Capacity (Water)
LP1111-B	Blue	200-225 Mpa	2.5	146.21	59-60cc
LP1111	White	200-225 Mpa	2.54	143.66	57-58cc

Composition		
Item Code	Cellulose Fiber Content (%)	Latex Content (%)
LP1111-B	100.0	0.0
LP1111	100.0	0.0

Contents			
Item Code	Total Size Inches	Quarter Folded Size, Inches (Avg.) ± 1/8"	Quantity
LP1111-B	12.5" x 14.4"	6 31/32 x 6 7/16	50 wipers/band 10 bands/case
LP1111	12.5" x 14.4"	6 31/32 x 6 7/16	50 wipers/band 10 bands/case

Physical Data is provided AS IS. Circumstances beyond Cancer Diagnostics' control may affect these values, including but not limited to, Storage Conditions, Ambient Humidity and Age.

Certified by: