



“The Cure for
the Common Lab”

Product Change Notification
Notification Date: May, 5 2016

PCN Number:

IA40001

Title:

Sprayer Cap Change for IA4000

Product Identification:

IA4000

Reason for Change:

Leak reduction during air shipment

Change Description:

The original sprayer cap has been replaced with a locking spray cap with an enhanced gasket. To lock the cap, turn the spray nozzle clock wise. To unlock the spray nozzle, turn the cap counter clock wise. The protective collar around the nozzle is removable but CDI recommends leaving it in place.

Identification Method to Distinguish Change:

New cap is semi translucent without a separate secondary cap.

Samples:

Available

Quantifiable Impact on Quality and Reliability:

None

Proposed First Ship Date:

5/5/2016