



Zero Xylene™

Odorless Xylene Substitute

#XSG1004 1 GAL (3.8L)

#XSI004 1 GAL (3.8L), 4/CS

DOT Proper Shipping Name:

3 / UN1268 / PGIII Petroleum Distillates, N.O.S.

(iso and cycloalkanes (C10-C11))

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____

Zero Xylene™ is a specially formulated clearant for tissue sections which does not contain xylenes or limonenes. Zero Xylene™ works well for routine H&E Staining and virtually all special stains and properly clears tissues sections following dehydration and facilitates penetration by the mounting medium for coverslipping.



DANGER

HAZARD STATEMENTS

Flammable liquid and vapor. May be fatal if swallowed and enters airways.

PRECAUTIONARY STATEMENTS - PREVENTION

Keep away from heat/sparks/open flames/hot surfaces. – No smoking. Keep container tightly closed. Ground / bond container and receiving equipment. Use explosion-proof electrical, ventilating, and lighting equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Avoid release to the environment. Wear protective gloves and eye / face protection.

PRECAUTIONARY STATEMENTS - RESPONSE

IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Do NOT induce vomiting.

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. If skin irritation occurs: Get medical advice/ attention.

IN CASE OF FIRE: Use water fog, foam, dry chemical or carbon dioxide (CO2) to extinguish.

STORAGE & DISPOSAL

Collect spillage. Store in a well-ventilated place. Keep cool. Store locked up.

Dispose of contents and container in accordance with local regulations.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

cancerdiagnostics.com

Made in USA