

# CASSETTE BLOCKS AND MICROSCOPE SLIDES STORAGE SYSTEMS



## Slide Genie™ Metal Slide System

Six drawer metal slide drawer holds up to 5,000 microscope slides. Stack 12 high on base unit for storage up to 60,000 slides. Each drawer has integrated separator which creates two rows of microscope slides per drawer.

SGM-6 Cabinet of 6 slide drawers  
SGM-B Base



**Secure Safety Latch**

## Microscope Slide Storage Box

Redesigned with enhanced sliding open/close locking mechanism. Foam cushion inserts, durable HDPE with hinged lid.

SB-25 25 slide capacity, white  
SB-50 50 slide capacity, white  
SB-100 100 slide capacity, white



**CANCER**  
DIAGNOSTICS, INC.

Cancer Diagnostics Inc.

Contact Information


Visit us online

4300 Emperor Blvd-400  
Durham, North Carolina 27703

Toll Free: (877) 846 5393  
Fax: (877) 817 1716

[cancerdiagnostics.com](http://cancerdiagnostics.com)

©2017 Cancer Diagnostics, Inc. All rights reserved. Printed in the USA



**EASY ORGANIZATION,  
CONVENIENT TRANSPORT & STORAGE**

# SuperLite™ Storage Units Microscope Slides



# SuperLite™ Storage Units Tissue Cassette



**SLIDE STORAGE ITEM# CPS-6**



PLASTIC INSERT

**SLIDE STORAGE ITEM# ACS-6**



CORRUGATED PARTITIONING

- 1 6 each, high impact polystyrene plastic trays, (7 1/2" wide x 13 1/4" deep x 2 5/8" high) /or Corrugated Partitioning
- 2 6 each, sturdy corrugated drawers with front labeling for easy retrieval
- 3 6 each, corrugated drawer enclosures for dust-proof protection
- 4 Reinforced Kraft Outer Carton for long term storage
- 5 Case dimensions: 27" wide x 16 1/2" deep x 10" high



\*Ships fully assembled

**CASSETTE STORAGE ITEM# CPC-6**



PLASTIC INSERT

**CASSETTE STORAGE ITEM# ACC-6**



CORRUGATED PARTITIONING

- 1 12 each, high impact polystyrene plastic trays, (7 1/2" wide x 13 1/4" deep x 1 3/4" high) /or Corrugated Partitioning
- 2 12 each, sturdy corrugated drawers with front labeling for easy retrieval
- 3 6 each, corrugated drawer enclosures for dust-proof protection
- 4 Reinforced Kraft Outer Carton for long term storage
- 5 Case dimensions: 27" wide x 16 1/2" deep x 10" high



\*Ships fully assembled

## Storage Files, Slide Filing System

CDI's corrugated interlocking storage files for microscope slides are engineered to be durable and functional. An outer shell acts as a stabilizer for added support. If desired, file boxes may also be used individually without being stacked. Designed to economically hold 2,000 slides per drawer cavity (10,000 slides per unit). Available flat, or fully assembled.

### Full/High Drawer:

- CG1111-X Fully Assembled CS/5 (1 Unit)
- CG1111-XF Flat/Unassembled CS/5 (1 Unit)
- CG1111-DBL-XF Flat/Unassembled CS/10 (2 Units)

### Original Half Drawer:

- CG1111 Fully Assembled CS/5 (1 Unit)
- CG1111-F Flat/Unassembled CS/5 (1 Unit)
- CG1111-DBL-F Flat/Unassembled CS/10 (2 Units)



## BlockFile™ Cassette Filing System

Lightweight and sturdy high impact poly construction, for long-term storage of tissue blocks. Six drawer cabinet accommodates up to 1,500 paraffin blocks.

CFS-1 Unit, BlockFile Cassette Filing System



## Storage Files, Cassette Filing System

CDI's corrugated interlocking files for tissue blocks are engineered to be durable and functional. An outer shell acts as a stabilizer to added support. If desired, file boxes may also be used individually without being stacked. Designed to hold 500 tissue blocks per drawer cavity (2,500 blocks per unit). Available flat, or fully assembled.

- CK1000 Fully Assembled CS/5 (1 Unit)
- CK1000-F Flat/Unassembled CS/5 (1 Unit)
- CK1000-DBL-F Flat/Unassembled CS/10 (2 Units)

