



Pap-Spray™

- Intended Use -

Pap-Spray™ is a polyethylene glycol formulation for rapid, high quality fixation of cytology cell spreads.

- General Information -

Pap-Spray™ is a specially formulated fixative that produces improved fixation of Pap smears and cell spreads. It combines coagulant fixation of ethanol with other enhanced fixation to produce better nuclear and cytoplasmic detail. Pap-Spray™ results in precise fixation of columnar and squamous cells while producing crisp, distinct cell membranes and nucleochromatin patterns.

- Application -

Simply spray a fine mist over your slide to preserve and protect the cells. A protective, thin film of water and alcohol-soluble polyethylene glycol forms over the sample, maximizing cell preservation and ensuring an optimum staining environment. This formulation is particularly applicable for microscopic examination of Pap smear.

- Packaging -

Convenient fine-mist sprayer. 5oz spray bottles for economical fixation. Also available in smaller 2oz bottles.

Catalog#	Volume
FX1086	5oz CS/12
FX1086-2	2oz CS/12

