



Acetone

#FX2101 1 GAL (3.8L)

#CX2101 5 GAL (19L) METAL PAIL

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____



DANGER
HAZARD STATEMENTS

Highly flammable liquid and vapor. Causes serious eye irritation.

PRECAUTIONARY STATEMENTS - PREVENTION

Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Wash body thoroughly after handling. Wear eye and face protection.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. If eye irritation persists: get medical attention.

IF ON SKIN: Take off immediately all contaminated clothing. Rinse skin with water/shower.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store in a well-ventilated place. Keep cool. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA