



Reagent Alcohol, 200 Proof

Specially Denatured Alcohol Formula 3-A. 95 parts
by volume and isopropyl 5 parts by volume

#AL0100 GAL (3.8L), CS/4

DOT Proper Shipping Name:

UN1170, Ethanol, 3, PG II

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____



DANGER

HAZARD STATEMENTS

Highly flammable liquid and vapor. Causes skin irritation. Causes eye irritation. May cause drowsiness or dizziness. Causes damage to organs

PRECAUTIONARY STATEMENTS - PREVENTION

Do not breathe mist/vapors/spray. Do not eat, drink or smoke when using this product. Ground/bond container and receiving equipment. Keep away from heat/sparks/open flames/hot surfaces-no smoking. Take precautionary measure against static discharge. Use only non-sparking tools. Use only outdoors or in a well-ventilated area. Wash thoroughly after handling. Wear protective gloves/eye protection/face protection.

PRECAUTIONARY STATEMENTS - RESPONSE

Call a poison center/doctor if you feel unwell. If exposed or concerned: Call a poison center/doctor. **IF IN EYES:** Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention **IF ON SKIN:** Take off immediately all

contaminated clothing. Rinse skin with water shower. Wash with plenty of water. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing

and wash it before reuse. **IF INHALED:** Remove person to fresh air and keep comfortable for breathing. **IN CASE OF FIRE:** Use consider water spray or fog, carbon dioxide, dry chemical powder, or alcohol resistant foam to extinguish.

STORAGE & DISPOSAL: Store in a well-ventilated place. Keep cool. Keep container tightly closed. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



CANCER
DIAGNOSTICS, INC.

4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA