



Michel's Solution

TRANSPORT MEDIA

#FX1046-500 500ML

#FX1046M 1000ML

#FX1046 1 GAL (3.8L)

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____

Michel's transport media is a solution that maintains isotonicity and pH of a tissue (7.0-7.2). It stabilizes proteins nicely for immunofluorescence.



WARNING **HAZARD STATEMENTS**

Harmful if swallowed.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Do not eat, drink, or smoke when using this product.

PRECAUTIONARY STATEMENTS - RESPONSE

IF SWALLOWED: Call a doctor if you feel unwell. Rinse mouth

STORAGE & DISPOSAL

Store in a well-ventilated place. Keep cool. Keep lid tightly closed. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com

Made in USA