



Date Received: _____ Date Opened: _____

A device that allows for accurate localization of suspect lesions in breast tissue (on radio-opaque grid). Grossing can also be performed directly on board. Non-Sterile. Single Use Only.


Insight-LLD™

Lesion Locating Device


Radiopaque Grid Board

REF #LLD16600

Box of 10



For In Vitro Diagnostic Use

59°F (15°C)  86°F (30°C)



4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com

MX2029 Insight-LLD Grid Board AW-2401395 Rev.B 2022-10-20