


Section 1. Identification

1.1 Product Identifier	Product Name: Propylene Glycol Product Number: SSC1114	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
		4300 Emperor Blvd. #400 Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 GHS Classification	Not classified as hazardous by the GHS.	
Classification: Directive 1999/45/EC	Not classified as hazardous by Directive 1999/45/EC. For full text of R- & S- phrases: see Section 15.	
2.2 Label Elements	Non hazardous mixture	
NFPA Scale: 0-4		
HMIS (U.S.A.) Scale: 0-4	HEALTH	0
	FLAMMABILITY	0
	PHYSICAL HAZARD	0
	PERSONAL PROTECTION	B
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.	

Section 3. Composition and Information on Ingredients

**3.2 Chemical Description:** Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Propylene Glycol	57-55-6	200-338-0	N/A	100

Section 4. First Aid Measures**4.1 Description of first aid measures****Eye Contact:** Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.**Skin Contact:** Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.**Inhalation:** If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.**Ingestion:** Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.**4.2 Most important symptoms and effects, both acute and delayed**

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	Hazardous combustion may release carbon monoxide, aldehydes.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures**6.1 Personal precautions, protective equipment and emergency procedures**

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage**7.1 Precautions for safe handling.**

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: N/A
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Clear
Odor	Odorless
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Strong oxidizing agent, strong acids, strong bases.
10.6 Hazardous Decomposition Materials	Oxides of carbon.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. May be irritating to skin and mucous membranes. Serious Eye Damage/Irritation. May be irritating to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods IATA: Not dangerous goods IMDG: Not dangerous goods ADR/RID: Not dangerous goods
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.

Section 16. Other Information
Revision Date: 2019-01-11


Notice to Reader:



To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Oil Red O Stain Product Number: SSC1100	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	4300 Emperor Blvd. #400 Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 GHS Classification	Carcinogenicity (Category 2) – H351	
Classification: Directive 1999/45/EC		Xn R40 EN: Harmful DA: Sundhedsskadelig DE: Gesundheitsschädlich EL: Επιβλαβές ES: Nocivo FR: Nocif IT: Nocivo NL: Schadelijk PT: Nocivo SV: Hälsoskadlig
	For full text of R- & S- phrases: see Section 15.	


2.2 Label Elements	Pictogram(s):  Signal word: Warning Hazard statement(s): H351 – Suspected of causing cancer Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes.										
NFPA Scale: 0-4											
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr><td>HEALTH</td><td>1</td></tr> <tr><td>FLAMMABILITY</td><td>0</td></tr> <tr><td>PHYSICAL HAZARD</td><td>0</td></tr> <tr><td>PERSONAL PROTECTION</td><td>B</td></tr> </table>	HEALTH	1	FLAMMABILITY	0	PHYSICAL HAZARD	0	PERSONAL PROTECTION	B		
HEALTH	1										
FLAMMABILITY	0										
PHYSICAL HAZARD	0										
PERSONAL PROTECTION	B										
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.										

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Propylene Glycol	57-55-6	200-338-0	N/A	≤ 100
Oil Red O	1320-06-5	215-295-3	 Warning 2, H351 Suspected of causing cancer	≤ 0.5

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious



person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	Hazardous combustion may release carbon monoxide, aldehydes.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: N/A
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection.

	Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Red
Odor	Odorless
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Strong oxidizing agents, strong acids, strong bases.
10.6 Hazardous Decomposition Materials	Oxides of carbon.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. May be irritating to skin and mucous membranes. Serious Eye Damage/Irritation. May be irritating to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed. Other: Benzidine-based dyes are recognized as possible carcinogens.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA, IMDG, ADR	DOT: Not dangerous goods ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods
14.2 UN Proper Shipping Name DOT, IATA, IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods

	DOT: Not dangerous goods
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods IATA: Not dangerous goods IMDG: Not dangerous goods ADR/RID: Not dangerous goods
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R40 EN: Limited evidence of a carcinogenic effect. DA: Mulighed for kræftfremkaldende effekt. DE: Verdacht auf krebserzeugende Wirkung. EL: Υποπτο καρκινογένεσης. ES: Posibles efectos cancerígenos. FR: Effet cancerogene suspecte ? preuves insuffisantes IT: Possibilità di effetti cancerogeni - prove insufficienti. NL: Carcinogene effecten zijn niet uitgesloten. PT: Possibilidade de efeitos cancerígenos. SV: Misstänks kunna ge cancer.
S phrases	S24/25 EN: Avoid contact with skin and eyes. DA: Undgå kontakt med huden og øjnene. DE: Berührung mit den Augen und der Haut vermeiden. EL: Αποφεύγετε επαφή με το δέρμα και με τα μάτια. ES: Evítese el contacto con los ojos y la piel. FR: Éviter le contact avec la peau et les yeux. IT: Evitare il contatto con gli occhi e con la pelle. NL: Aanraking met de ogen en de huid vermijden. PT: Evitar o contacto com a pele e os olhos. SV: Undvik kontakt med huden och ögonen. S36/37/39 EN: Wear suitable protective clothing, gloves and eye/face protection. DA: Brug særligt arbejdstøj, egnede beskyttelseshandsker og -briller/ansigtsskærm. DE: Bei der Arbeit geeignete Schutzkleidung, Schutzhandschuhe und Schutzbrille/Gesichtsschutz tragen. EL: Φοράτε κατάλληλη προστατευτική ενδυμασία, γάντια και συσκευή προστασίας ματιών/προσώπου ES: Usense indumentaria y guantes adecuados y protección para los ojos/la cara.

	<p>FR: Porter un vêtement de protection approprié, des gants et un appareil de protection des yeux/du visage.</p> <p>IT: Usare indumenti protettivi e guanti adatti e proteggersi gli occhi/la faccia.</p> <p>NL: Draag geschikte beschermende kleding, handschoenen en een beschermingsmiddel voor de ogen/het gezicht.</p> <p>PT: Usar vestuário de protecção, luvas e equipamento protector para os olhos /face adequados.</p> <p>SV: Använd lämpliga skyddskläder, skyddshandskar samt skyddsglasögon eller ansiktsskydd.</p>
--	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Section 16. Other Information

Revision Date: 2019-01-11


Notice to Reader:



To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Hematoxylin, Mayer's (Lillie's Modification) Product Number: SSC1059	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	4300 Emperor Blvd. #400 Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 GHS classification	Acute toxicity, Oral (Category 4) – H302 Skin irritation (Category 2) – H315 Eye irritation (Category 2A) – H319
Classification: Directive 1999/45/EC	 <p>Xi, R36/38; S2425, S26, S45 EN: Irritant DA: Lokalirriterende DE: Reizend EL: Ερεθιστικό ES: Irritante FR: Irritant IT: Irritante NL: Irriterend PT: Irritante SV: Irriterande</p> <p>For full text of R- & S- phrases: see Section 15.</p>





2.2 Label Elements	<p>Pictogram(s): </p> <p>Signal word: Warning</p> <p>Hazard statement(s): H302 – Harmful if swallowed. H315 – Causes skin irritation H320 – Causes serious eye irritation</p> <p>Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes. P301 + P312 - IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.</p>										
NFPA Scale: 0-4											
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr> <td>HEALTH</td> <td>1</td> </tr> <tr> <td>FLAMMABILITY</td> <td>0</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>B</td> </tr> </table>	HEALTH	1	FLAMMABILITY	0	PHYSICAL HAZARD	0	PERSONAL PROTECTION	B		
HEALTH	1										
FLAMMABILITY	0										
PHYSICAL HAZARD	0										
PERSONAL PROTECTION	B										
2.3 Other Hazards	<p>PBT: This mixture does not contain any substances that are assessed to be a PBT.</p> <p>vPvB: This mixture does not contain any substances that are assessed to be a vPvB.</p>										



Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Glycerin	56-81-5	200-289-5	n/a	≤ 30
Aluminum Sulfate	7784-31-8	n/a	 Warning. 2 H315 Causes skin irritation. 3 H335 May cause respiratory irritation.  Danger. 2 H318 Causes serious eye damage.  Warning. 2 H411 Toxic to aquatic life with long lasting effects.	≤ 5
Acetic Acid	64-19-7	200-580-7	 Warning. 3, H226	≤ 2

			 Flammable liquid and vapour. Danger. 1A H314 Causes severe skin burns and eye damage.	
Hematoxylin	517-28-2	208-237-3	 Warning. 4 H302 Harmful if swallowed. 2 H315 Causes skin irritation. 2 H319 Causes serious eye irritation. 3 H335 May cause respiratory irritation.	≤ 1

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.
 Wash thoroughly after handling.
 Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.
 Keep container tightly closed.
 Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Acetic acid, glacial; ACGIH: 15ppm STEL (37 mg/m ³ STEL) 10 ppm TWA (25 mg/m ³ TWA); NIOSH: 15 ppm STEL (37 mg/m ³ STEL) 10 ppm TWA (25 mg/m ³ TWA); OSHA: 10 ppm TWA (25 mg/m ³ TWA)
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Purple
Odor	Odorless
Odor Threshold	Unknown
pH	2.4 ± 0.3
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown

Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Reactive with strong oxidizing agents, reducing agents, metals, acids, and alkalis.
10.6 Hazardous Decomposition Materials	Carbon Monoxide. Carbon Dioxide. May release irritating or toxic fumes under fire conditions.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	<p>Acute Toxicity. No relevant data available</p> <p>Skin Corrosion/Irritation. Irritating to skin and mucous membranes.</p> <p>Serious Eye Damage/Irritation. Irritating to eyes.</p> <p>Respiratory or skin sensitization. No relevant data available.</p> <p>Germ Cell Mutagenicity. No relevant data available.</p>
Carcinogenicity.	<p>International Agency for Research on Cancer (IARC). None of the components are listed.</p> <p>National Toxicology Program (NTP). None of the components are listed.</p>

Section 12. Ecological Information

12.1 Toxicity	<p>Fish: No relevant studies identified.</p> <p>Crustacea: No relevant studies identified.</p> <p>Algae/Aquatic Plants: No relevant studies identified.</p> <p>Other Organisms: Acetic Acid: LD50 Oral, Rat, 3310 mg/kg. LD50 Dermal, Rabbit, 1060 mg/kg.</p>
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.

12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods. ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods.
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods. IATA: Not dangerous goods. IMDG: Not dangerous goods. ADR/RID: Not dangerous goods.
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.
14.7 Transport in Bulk According to Annex II of MARPOL73/78 and the IBC Code.	This product is not intended to be transported in bulk as defined in Annex II of MARPOL73/78 and the IBC Code.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous	None of the components in this mixture are listed.

Substances; Section 355	
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	This mixture does not contain a chemical known to the state of California to cause cancer.
Right To Know Components	<p>Acetic Acid: California Massachusetts Minnesota New Jersey Pennsylvania Rhode Island</p> <p>Glycerin: Illinois Massachusetts Minnesota Pennsylvania Rhode Island</p> <p>Acetic Acid is listed on Canada's Ingredient Disclosure List.</p>
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	<p>R36/38 EN: Irritating to eyes and skin. DA: Irriterer øjnene og huden DE: Reizt die Augen und die Haut. EL: Ερεθίζει τα μάτια και το δέρμα. ES: Irrita los ojos y la piel. FR: Irritant pour les yeux et la peau. IT: Irritante per gli occhi e la pelle. NL: Irriterend voor de ogen en de huid. PT: Irritante para os olhos e pele. SV: Irriterar ögonen och huden.</p>
S phrases	<p>S24/25 EN: Avoid contact with skin and eyes. DA: Undgå kontakt med huden og øjnene. DE: Berührung mit den Augen und der Haut vermeiden. EL: Αποφεύγετε επαφή με το δέρμα και με τα μάτια. ES: Evítese el contacto con los ojos y la piel. FR: Éviter le contact avec la peau et les yeux. IT: Evitare il contatto con gli occhi e con la pelle. NL: Aanraking met de ogen en de huid vermijden. PT: Evitar o contacto com a pele e os olhos. SV: Undvik kontakt med huden och ögonen.</p> <p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen.</p>

	<p>PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista.</p> <p>SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S45</p> <p>EN: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).</p> <p>DA: Ved ulykkestilfælde eller ved ildebefindende er omgående lægebehandling nødvendig; vis etiketten, hvis det er muligt.</p> <p>DE: Bei Unfall oder Unwohlsein sofort Arzt zuziehen (wenn möglich, dieses Etikett vorzeigen).</p> <p>EL: Σε περίπτωση ατυχήματος ή αν αισθανθείτε αδιαθεσία ζητήστε άμεσα ιατρική συμβουλή (δείξτε την ετικέτα αν είναι δυνατό).</p> <p>ES: En caso de accidente o malestar, acúdase inmediatamente al médico (si es posible, muéstresele la etiqueta).</p> <p>FR: En cas d'accident ou de malaise, consulter immédiatement un medecin (si possible lui montrer l'etiquette).</p> <p>IT: In caso di incidente o di malessere consultare immediatamente il medico (se possibile, mostrargli l'etichetta).</p> <p>NL: Bij een ongeval of indien men zich onwel voelt, onmiddellijk een arts raadplegen (indien mogelijk hem dit etiket tonen).</p> <p>PT: Em caso de acidente ou de indisposição, consultar imediatamente o médico (se possível mostrar-lhe o rótulo).</p> <p>SV: Vid olycksfall, illamående eller annan påverkan, kontakta omedelbart läkare. Visa om möjligt etiketten.</p>
--	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Section 16. Other Information

Revision Date: 2019-01-11

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.