



CANCER
DIAGNOSTICS, INC.

PureAffin[®] 1

Tissue Infiltration Medium

REF PARA1

1 kg x 8

A purified paraffin specially formulated to provide ideal infiltration with all types of tissue specimens. Tissue Infiltration Medium is compatible with all industry paraffins.

Use with open or closed tissue processors.

A Compound of Purified Paraffin and additives.

Stable at elevated temperatures.

For use as an infiltration and embedding medium in histology.

Melting Temp: 56-57°C (133-135F)

WARNING!

May cause respiratory irritation, May cause damage to organs through prolonged or repeated exposure. Do not breathe dust/fume/gas/mist/vapors/spray. Get medical attention/advice if you feel unwell. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Store in a well-ventilated place. Keep container tightly closed, Dispose of contents/container to an approved waste disposal plant.

COMPOSITION/INFORMATION ON INGREDIENTS :

Paraffin Waxes and Hydrocarbon Waxes CAS# 8002-74-2, Polymer Wax CAS# 68527-08-2

Case reorder # PARA1

For In Vitro Diagnostic Use



59°F (15°C) 86°F (30°C)

EC REP

mil Europa GmbH Langenhagener Str. 71 D-30855 Langenhagen Germany +49-511-3908 9530



Cancer Diagnostics, Inc. 4300 Emperor Blvd - 400 Durham, NC 27703 USA

MX1 124 PureAffin 1, INFILTRATION PARAFFIN, BK/GCS
AW-PARA1 Rev.A 2015-02-20

IVD CE MADE IN USA

Order Toll Free: (877) 846 5393

Fax: (877) 817 1716

www.cancerdiagnostics.com