



CANCER
DIAGNOSTICS, INC.

PureAffin[®]

Tissue Embedding/Infiltration Medium
Sectioning with minimal compression

REF PARA2

1 kg x 8

Use for both processing and embedding or with PureAffin 1 Infiltration Medium as part of two step system. Compatible with all industry paraffins. Blended to provide thorough impregnation of tissues and facilitate excellent sectioning characteristics.

Use with open or closed tissue processors.
A Compound of Purified Paraffin and additives.
Stable at elevated temperatures.

For use as an infiltration and embedding medium in histology.

Melting Temp: 56-57°C (133-135F)

WARNING!

May cause respiratory irritation, May cause damage to organs through prolonged or repeated exposure. Do not breathe dust/fume/gas/mist/vapors/spray. Get medical attention/advice if you feel unwell. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Store in a well-ventilated place. Keep container tightly closed. Dispose of contents/container to an approved waste disposal plant

COMPOSITION/INFORMATION ON INGREDIENTS

Paraffin Waxes and Hydrocarbon Waxes CAS# 8002-74-2, Polymer Wax CAS# 68527-08-2

Case reorder # PARA2

For In Vitro Diagnostic Use



EC REP mdi Europa GmbH Langenhagener Str. 71 D-30855 Langenhagen Germany +49-511-3908 9530

Cancer Diagnostics, Inc. 4300 Emperor Blvd - 400 Durham, NC 27703 USA

MEX1125 PureAffin 2, INFILTRATION AND EMBEDDING PARAFFIN, 8KGCS IVD CE MADE IN USA
AW-PARA2 Rev. A 2015-02-20

Order Toll Free: (877) 846 5393

Fax: (877) 817 1716

www.cancerdiagnostics.com