



PureAffin[®]

Ultra-Pure Infiltration Paraffin

REF EEPAR1

1 kg x 8

Formulated with DMSO for enhanced infiltration. Use for both infiltration and embedding, or as an excellent infiltration paraffin in your two-step paraffin system. Designed to infiltrate fatty tissue with ease. Lends a degree of elasticity to the tissue.

For use as an infiltration and embedding medium in histology.

Melting Temp: 56-57 C (132.8 – 134.6 F)

WARNING!

May cause respiratory irritation. May cause damage to organs through prolonged or repeated exposure. Do not breathe dust/fumes/gas/mist/vapors/spray. Get medical attention/advice if you feel unwell. IF INHALED: Remove victim to fresh air and keep them at rest in a position comfortable for breathing. Store in a well-ventilated place. Keep container tightly closed. Dispose of contents/container to an approved waste disposal plant.

COMPOSITION/INFORMATION ON INGREDIENTS:

Paraffin Waxes and Hydrocarbon Waxes CAS# 8002-74-2, Polymer Wax CAS# 68527-08-2

Case reorder # EEPAR1

LOT



For In Vitro Diagnostic Use



59°F (15°C)  86°F (30°C)

EC REP mdi Europa GmbH Langenhagener Str.71 D-30855 Langenhagen Germany +49-51 1-3908 9530

 Cancer Diagnostics, Inc. 4300 Emperor Blvd • 400 Durham, NC 27703 USA

MX1131 PureAffin I INFILTRATION PARAFFIN, 8KG/CS AW-EEPAR1 Rev.A 2016-4-15
AW-EEPAR1 Rev.A 2016-4-15

IVD CE MADE IN USA

Order Toll Free: (877) 846 5393

Fax: (877) 817 1716

cancerdiagnostics.com