Complies with EC no. 1907/2006

Date of Issue: 10/27/1999

Date of Revision: 09/19/2018

# Safety Data Sheet (SDS)

## **Section 1: Chemical Product and Company Identification**

Cat#: 1601

Part Name: CitriSolv

Supplier: Decon Laboratories Inc.

460 Glennie Circle King of Prussia, Pa 19406

SDS Telephone # (610) 755-0800

Identified uses: Laboratory use

Email Contact: cveloski@deconlabs.com

### **Emergency Telephone Numbers**

US Chemtrec: (800) 424-9300 Canada: (703) 527-3887

#### **Section 2: Hazards Identification:**

#### **Hazard Overview**

Classification of the substance or mixture

## Signal Word: WARNING









### Hazard and Precautionary Statements

H226 Flammable liquid and vapor H351: Suspected of causing cancer H315 Causes skin irritation.

H317: May cause an allergic skin reaction

H400: Very toxic to aquatic life

H410: Very toxic to aquatic life with long lasting effects

#### **Precautionary statements**

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

Use explosion-proof electrical/ventilating/lighting/equipment.

Avoid breathing dust/fume/gas/mist/vapors/spray

Wear protective gloves / eye protection / face protection. Wear protective gloves.

Ground/bond container and receiving equipment.

Keep container tightly closed.

Use only non-sparking tools.

Avoid release to the environment.

Take precautionary measures against static discharge.

Wash thoroughly after handling.

Contaminated work clothing must not be allowed out of the workplace.

Obtain special instructions before use.

Do not handle until all safety precautions have been read and understood.

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/ shower.

Specific treatment (see on this label).

Wash contaminated clothing before reuse.

IF exposed or concerned: Get medical advice/attention.

If skin irritation or rash occurs: Get medical advice/attention.

In case of fire: Use for extinction: CO2, powder or water spray.

Collect spillage.

Take off contaminated clothing and wash it before reuse.

Store locked up.

Store in a well-ventilated place. Keep cool.

Dispose of contents/container in accordance with local/regional/national/international regulations.

# Section 3: Composition/Information on ingredients

Common Name	Chemical Name	CAS	EINECS	Weight %
		Number	No.	
d-Limonene	(R)-1-methyl-4-(1- methylenthenyl)- Cyclohexene	5989-27-5	277-813- 5	> 98.5

#### Section 4: First Aid Measures

- · Description of first aid measures
- After inhalation:

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

- · After skin contact: Immediately wash with water and soap and rinse thoroughly.
- · After eye contact: Rinse opened eye for several minutes under running water.
- · After swallowing: If symptoms persist consult doctor.
- · Information for doctor:
- Most important symptoms and effects, both acute and delayed No further relevant information available.
- Indication of any immediate medical attention and special treatment needed No further relevant information available.

### **Section 5: Fire-Fighting Measures**

- · Extinguishing media
- · Suitable extinguishing agents: CO2, sand, extinguishing powder. Do not use water.
- · For safety reasons unsuitable extinguishing agents: Water with full jet
- · Special hazards arising from the substance or mixture No further relevant information available.
- Advice for firefighters

Protective equipment: No special measures required.

#### Section 6: Accidental Release measures

Any information listed below is to be considered in addition to internal guidelines for isolation of spill, containment of spill, removal of ignition sources from immediate area, and collection for disposal of spill by trained, properly protected clean up personnel.

- · Personal precautions, protective equipment and emergency procedures
- Wear protective equipment. Keep unprotected persons away.
- · Environmental precautions:

Do not allow product to reach sewage system or any water course. Inform respective authorities in case of seepage into water course or sewage system. Do not allow to enter sewers/ surface or ground water.

· Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust). Ensure adequate ventilation.

Do not flush with water or aqueous cleansing agents

· Reference to other sections

See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment. See Section 13 for disposal information.

## **Section 7: Handling and Storage**

#### Handling:

#### Precautions for safe handling

Ensure good ventilation/exhaustion at the workplace.

Prevent formation of aerosols.

Information about protection against explosions and fires:

Keep ignition sources away – Do not smoke.

Protect against electrostatic charges.

Conditions for safe storage, including any incompatibilities

Storage:

Requirements to be met by storerooms and receptacles: No special requirements.

Information about storage in one common storage facility: Not required Further information about storage conditions: Keep receptacle tightly sealed.

**Specific end use(s)** No further relevant information available.

# **Section 8: Exposure Controls/ Personal Protection**

- · Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
- · Components with limit values that require monitoring at the workplace: Not required.
- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls

#### · Personal protective equipment:

#### General protective and hygienic measures:

Keep away from foodstuffs, beverages and feed. Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the skin.

Avoid contact with the eyes and skin.

#### · Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· Protection of hands:



#### Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

#### Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

#### · Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection:



#### Tightly sealed goggles

Chemical-resistant gloves should be worn whenever this material is handled. The glove material has to be impermeable and resistant to the product. Gloves should be removed and replaced immediately if there is any indication of degradation or chemical breakthrough. Rinse and remove gloves immediately after use. Wash hands with soap and water. All glove recommendations presume that the risk of exposure is through splash and not intentional immersion of the hands into the product. Since glove permeation data does not exist for this material, no recommendation for the glove material can be given for the product. Permiation data must be obtained from the glove manufacturer to determine if the glove is suitable for the task.

### **Section 9: Physical and Chemical Properties**

Information on books abusingles and a	homical proportion	
<ul> <li>Information on basic physical and c</li> <li>General Information</li> </ul>	enemical properties	
· Appearance:		
Form:	Liquid	
Color:	Colorless	
· Odor:	Characteristic	
· Odour threshold:	Not determined.	
· pH-value:	Not determined.	
· Change in condition		
Melting point/Melting range:	-74.3 °C (-102 °F)	
Boiling point/Boiling range:	ab173 °C (ab343 °F)	
· Flash point:	43 °C (109 °F)	
· Flammability (solid, gaseous):	Not applicable.	
· Ignition temperature:	255 °C (491 °F)	
· Decomposition temperature:	Not determined.	
· Auto igniting:	Not determined.	
<ul><li>Danger of explosion: air/</li></ul>	Product is not explosive. However, formation of explosive vapor mixtures are possible.	
· Explosion limits:		
Lower:	0.7 Vol %	
Upper:	6.1 Vol %	
Vapor pressure at 20 °C (68 °F):	2.3 hPa (2 mm Hg)	
· Density at 20 °C (68 °F):	0.848 g/cm³ (7.077 lbs/gal)	
· Relative density	Not determined.	
· Vapour density	Not determined.	
· Evaporation rate	Not determined.	
· Solubility in / Miscibility with		
Water:	Not miscible or difficult to mix.	
· Partition coefficient (n-octanol/wate	er): Not determined.	
· Viscosity:		
Dynamic:	Not determined.	
Kinematic:	Not determined.	
Organic solvents:	99.8 %	
VOC content:	99.8 %	
Other to farme at the	846.0 g/l / 7.06 lb/gl	
· Other information	No further relevant information available.	

# **Section 10: Stability and Reactivity:**

Reactivity
Chemical stability

#### · Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

## **Section 11: Toxicological Information**

- · Information on toxicological effects
- · Acute toxicity:

LD/LC50 values that are relevant for classification:			
5989-27-5 (R)-p-mentha-1,8-diene			
Oral	LD50	4400 mg/kg (rat)	

#### · Primary irritant effect:

- · on the skin: Irritant to skin and mucous membranes.
- · on the eye: No irritating effect.
- · Sensitization: Sensitization possible through skin contact.
- · Additional toxicological information:
- · Carcinogenic categories

IARC (International Agency for Research on Cancer)				
5989-27-5	(R)-p-mentha-1,8-diene	3		
25013-16-5	tert-butyl-4-methoxyphenol	2B		
NTP (National Toxicology Program)				
25013-16-5	tert-butyl-4-methoxyphenol	R		
OSHA-Ca (Occupational Safety & Health Administration)				
Substance is not	listed.			

## **Section 12: Ecological Information**

- Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Ecotoxical effects:
- · Remark: Very toxic for fish
- · Additional ecological information:
- · General notes:

Water hazard class 2 (Assessment by list): hazardous for water

Do not allow product to reach ground water, water course or sewage system.

Danger to drinking water if even small quantities leak into the ground.

Also poisonous for fish and plankton in water bodies.

Very toxic for aquatic organisms

- · Results of PBT and vPvB assessment
- · PBT: Not applicable.

· vPvB: Not applicable.

Other adverse effects No further relevant information available.

### **Section 13: Disposal Considerations**

- · Waste treatment methods
- · Recommendation:

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

- · Uncleaned packagings:
- · Recommendation: Disposal must be made according to official regulations.

### **Section 14: Transportation Information**

Proper Shipping Name:

Cleaning Compound Liquid (Non-Hazardous)

### **Section 15: Regulatory Information**

- · Safety, health and environmental regulations/legislation specific for the substance or mixture
- · Sara
  - · Section 355 (extremely hazardous substances):

Substance is not listed.

· Section 313 (Specific toxic chemical listings):

Substance is not listed.

· TSCA (Toxic Substances Control Act):

Substance is listed.

· Proposition 65

This product has been classified in accordance with the hazard criteria of the CPR (Controlled Products Regulations) and this SDS (Safety Date Sheet) contains all information required by CPR.

#### Section 16: Other Information

Date of Issue: 10/27/1999 Date of Revision: 09/19/2018

Decon Laboratories, Inc. provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. Individuals receiving this information must exercise their independent judgment in determining its appropriateness for a particular purpose. Decon Laboratories, Inc. makes no representations or warranties, either expressed or implied of merchantability, fitness for particular purposes with respect to the information set forth herein or to which the information refers. Accordingly, Decon Laboratories, Inc. will not be responsible for damages resulting from the use of or reliance upon this information.

**End of Safety Data Sheet**