



# Reagent Alcohol, 95%

Specially Denatured Alcohol Formula 3-A.  
diluted with 5 parts deionized water  
by volume, 190 Proof

#AL5095 5 GAL (3.8L) CUBE

DOT Proper Shipping Name:

**UN1170 Ethanol Solutions, 3, PG II**

**For In Vitro Diagnostics Use**

Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_



## DANGER

### HAZARD STATEMENTS

Highly flammable liquid and vapor. Causes skin irritation. Causes eye irritation. May cause drowsiness or dizziness. Causes damage to organs

### PRECAUTIONARY STATEMENTS - PREVENTION

Do not breathe mist/vapors/spray. Do not eat, drink or smoke when using this product. Ground/bond container and receiving equipment. Keep away from heat/sparks/open flames/hot surfaces-no smoking. Take precautionary measure against static discharge. Use only non-sparking tools. Use only outdoors or in a well-ventilated area. Wash thoroughly after handling. Wear protective gloves/eye protection/face protection.

### PRECAUTIONARY STATEMENTS - RESPONSE

Call a poison center/doctor if you feel unwell. If exposed or concerned: Call a poison center/doctor. **IF IN EYES:** Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention **IF ON SKIN:** Take off immediately all contaminated clothing. Rinse skin with water shower. Wash with plenty of water. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. **IF INHALED:** Remove person to fresh air and keep comfortable for breathing. **IN CASE OF FIRE:** Use consider water spray or fog, carbon dioxide, dry chemical powder, or alcohol resistant foam to extinguish.

**STORAGE & DISPOSAL:** Store in a well-ventilated place. Keep cool. Keep container tightly closed. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

**NOTICE:** Read Safety Data Sheet for additional information.

**EMERGENCY CONTACT: Chemtrec 1-800-424-9300**



CANCER  
DIAGNOSTICS, INC.

4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

**Made in USA**