



Clear Away Citrus

REF #FX2180 1 Gallon (3.8L)

LOT



For In Vitro Diagnostics Use

46°F (8°C)
35°F (2°C)

Date Received: _____ Date Opened: _____

Clear Away Citrus is a citrus clearing agent formulated specifically to replace xylene and toluene in histology and cytology laboratory procedures. It is less oily than other citrus solvent substitutes, dries quickly, and does not leave a residue. Clear Away Citrus penetrates tissue rapidly without over-hardening and is safe to use with all industry tissue processors. Clear Away Citrus is biodegradable and offers the desirable characteristics of a clearing solvent without the harmful effects of aromatic hydrocarbons such as xylene and other toluene-based products.



DANGER HAZARD STATEMENTS

Causes skin and serious eye irritation. May cause an allergic skin reaction. Flammable liquid and vapor.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Avoid breathing dust, vapors. Contaminated work clothing must not be allowed out of the workplace. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. If eye irritation persists: get medical attention.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off immediately all contaminated clothing and wash it before reuse.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store in a well-ventilated place at room temperature. Keep cool. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**

MX1646 ClearAway Citrus 1 Gallon AW-6073233 Rev.B 2021-03-12



Clear Away Citrus

REF #FX2181 | 1 Gallon (3.8L) CS/4

LOT



For In Vitro Diagnostics Use

46°F (8°C)
35°F (2°C)

Date Received: _____ Date Opened: _____

Clear Away Citrus is a citrus clearing agent formulated specifically to replace xylene and toluene in histology and cytology laboratory procedures. It is less oily than other citrus solvent substitutes, dries quickly, and does not leave a residue. Clear Away Citrus penetrates tissue rapidly without over-hardening and is safe to use with all industry tissue processors. Clear Away Citrus is biodegradable and offers the desirable characteristics of a clearing solvent without the harmful effects of aromatic hydrocarbons such as xylene and other toluene-based products.



DANGER HAZARD STATEMENTS

Causes skin and serious eye irritation. May cause an allergic skin reaction. Flammable liquid and vapor.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Avoid breathing dust, vapors. Contaminated work clothing must not be allowed out of the workplace. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. If eye irritation persists: get medical attention.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off immediately all contaminated clothing and wash it before reuse.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store in a well-ventilated place at room temperature. Keep cool. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**

MX1647 ClearAway Citrus Gallons, 4 AW-6073234 Rev.B 2021-03-12