



# Cytology Clarifier

#FX2102 1000ML

For In Vitro Diagnostics Use

Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_

Cytology Clarifier™ is formulated to differentiate hematoxylin in regressive cytology H&E staining procedures. Nuclei demonstrate precisely stained nuclear chromatin, nucleoli are discrete and nuclear membranes are distinct. Non-nuclear structures show excellent staining with cytology counterstains, with ideal variations of staining of erythrocytes and nucleated cell cytoplasm.

Mucin is not stained with hematoxylin following differentiation with Cytology Clarifier™. This is the optimum differentiator for use with CDI's Harris Hematoxylin.



## **WARNING** **HAZARD STATEMENTS**

Causes skin irritation. Causes eye irritation.

## **PRECAUTIONARY STATEMENTS - PREVENTION**

Wash body thoroughly after handling. Wear protective gloves.

## **PRECAUTIONARY STATEMENTS - RESPONSE**

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. If eye irritation persists: get medical attention.

**IF ON SKIN:** Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

## **STORAGE & DISPOSAL**

Store in a well-ventilated, cool place. Keep lid tightly closed. Dispose of contents/container to an approved waste disposal plant.

**NOTICE:** Read Safety Data Sheet for additional information.

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)

**Made in USA**