



# Thin Layer

## Non-Gyn

### Cell Preserve Solution, 20mL

REF #1NGYN-100V CS/100

LOT



FDX1268 Thin Layer Cell Preserve Solution Val AW-14104 Rev.B 2020-04-13

### For In Vitro Diagnostics Use

Preservative solution for cell support during transport and preparation of slides. For use with liquid based cytology method. For Non-Gyn cytological specimens only.



#### DANGER - HAZARD STATEMENTS

Causes skin irritation. Causes serious eye irritation. Causes damage to organs. May cause damage to organs through prolonged or repeated exposure. Flammable liquid and vapor. Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Do not breathe vapors. Do not eat, drink or smoke when using this product. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Full Labeling information for Thin Layer Cell Preserve Solution located on outer packaging. Store in product packaging.

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



CANCER  
DIAGNOSTICS INC.

59°F (15°C)



86°F (30°C)

4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)