



Nail-RX™

KERATIN SOFTENING AGENT

REF #NL0500 500ML

For In Vitro Diagnostics Use



Date Received: _____ Date Opened: _____

Use pre-processing to soften nail clippings. Nail-RX™ nail prep solution softens brittle nails, making them easier to cut and process, while helping to keep the specimen from washing off the slide. Soak for 1 hour and process as usual. (Longer depending on size of specimen.)



DANGER. CONTAINS FORMALDEHYDE HAZARD STATEMENTS

Harmful if inhaled. Causes severe skin burns and serious eye damage. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction, genetic defects, cancer. May cause damage to organs (lungs, nose) through single, prolonged, or repeated exposure. May be corrosive to metals.

PRECAUTIONARY STATEMENTS - PREVENTION

Avoid breathing dust, vapors. Use only outdoors or in a well-ventilated area. Do not breathe dusts, mists or vapors. Wash body thoroughly after handling. Wear protective gloves and clothing. Wear eye and face protection. Wear NIOSH approved respiratory protection. Contaminated work clothing must not be allowed out of the workplace. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Do not eat, drink or smoke when using this product. Keep only in original container.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing.

IF ON SKIN: Rinse skin with water/shower. Take off immediately all contaminated clothing and wash it before reuse. Immediately call a doctor. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. If exposed or concerned: get medical attention/call a doctor.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF INHALED: If breathing is difficult, remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: Call a doctor. Call a doctor if you feel unwell.

STORAGE & DISPOSAL

Store locked up. Store in corrosive resistant container with a resistant inner liner. Absorb spillage to prevent material damage. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**

MS2271 104-RX 01/14 0118 Rev. B 2020/04/29