



**PUREVIEW®**  
Hematoxylin and Eosin  
Staining System

**FOCUS-AQ™**

(Aqueous) Concentrate

**REF**

16703-C 500mL each

16703-C4 500mL, 4/case



Order Toll Free: (877) 846 5393

Fax: (877) 817 1716

[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)

 Cancer Diagnostics, Inc.  
4300 Emperor Blvd-400.  
Durham, NC 27703

NET 500 ML

Made in USA

**DIRECTIONS:** Focus AQ is a concentrated reagent and should be diluted before use. To obtain a working solution, pour two measures of concentrate into 500mL of distilled water (30mL). One measure equals fluid filled to first line in bottle cap. A small margin of error produced from using the cap for measurement is insignificant to the performance, pH, and strength of the reagent.

For in vitro diagnostic use only.  
Store at room temperature (20°-30°C).

**DANGER**

**HAZARD STATEMENT(S):** Causes skin irritation. Causes eye irritation

**PRECAUTIONARY STATEMENT(S):** Wash body thoroughly after handling. Wear protective gloves.

**IF ON SKIN (OR HAIR):** Wash with plenty of water. Specific treatment, see first aid section on the SDS. If skin irritation or rash occurs: Get medical attention. Take off all contaminated clothing and wash it before reuse.

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical attention.

**NOTICE:** Read Safety Data Sheet (SDS) for additional information.

24 Hour Emergency Contact:  
CHEMTREC 800-424-9300  
USA and Canada



Date Received: \_\_\_\_\_

Date Opened: \_\_\_\_\_