



Pureview® Focus-AQ™

Working Solution

REF #16703 | 1 Gallon (3.8L)

LOT



For In Vitro Diagnostics Use



Date Received: _____ Date Opened: _____

Pureview® Focus® is a clarifier that removes non-specific staining and background haze. Focus is formulated to sharpen nuclear detail producing a finishing effect that optimizes nuclear detail.



WARNING **HAZARD STATEMENTS**

Causes skin irritation. Causes serious eye irritation

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing If eye irritation persists: Get medical attention.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

STORAGE & DISPOSAL

Store at room temperature. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

MX1533 PureView Focus-AQ RTU 1 Gallon AWW-5672864 Rev.B 2020-10-16



Pureview® Focus-AQ™

Working Solution

REF #16703-4 | 1 Gallon (3.8L) CS/4

LOT



For In Vitro Diagnostics Use



Date Received: _____ Date Opened: _____

Pureview® Focus® is a clarifier that removes non-specific staining and background haze. Focus is formulated to sharpen nuclear detail producing a finishing effect that optimizes nuclear detail.



WARNING **HAZARD STATEMENTS**

Causes skin irritation. Causes serious eye irritation

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing If eye irritation persists: Get medical attention.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

STORAGE & DISPOSAL

Store at room temperature. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

MX1534 PureView Focus-AQ RTU Gallons, 4 AW-5672865 Rev.B 2020-10-16