



PUREVIEW[®]
Hematoxylin and Eosin
Staining System

Eosin 615[™]

with Phloxine

REF

16601-PH 500mL each
16601-PH4 500mL, 4/case



Order Toll Free: (877) 846 5393

Fax: (877) 817 1716

www.cancerdiagnostics.com

 Cancer Diagnostics, Inc.
4300 Emperor Blvd-400.
Durham, NC 27703

NET 500 ML
Made in USA

For in vitro diagnostic use only.
Store at room temperature (20°-30°C).

DANGER

HAZARD STATEMENT(S): Causes skin irritation. Causes serious eye irritation. Causes damage to organs (respiratory system). May cause damage to organs (respiratory system, central nervous system, liver, blood) through prolonged or repeated exposure. Highly flammable liquid and vapor.

PRECAUTIONARY STATEMENT(S): Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Do not breathe dust, vapors. Do not eat, drink, or smoke when using this product. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

IF ON SKIN (OR HAIR): Wash with plenty of water. Specific treatment, see first aid section on the SDS. If skin irritation or rash occurs: Get medical attention. Take off all contaminated clothing and wash it before reuse.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical attention. If exposed or concerned: call a doctor. Call a doctor if you feel unwell.

IN CASE OF FIRE: Use water, dry chemical, CO₂ or foam to extinguish.

STORAGE: Store locked up. Store in a well-ventilated place. Keep cool.

DISPOSAL: Dispose of contents/container in accordance with local regulations.

NOTICE: Read Safety Data Sheet (SDS) for additional information.

24 Hour Emergency Contact:
CHEMTREC 800-424-9300
USA and Canada



Date Received: _____

Date Opened: _____