



10% Alcoholic Formalin

Fixative

REF #FX1027 1 Gallon (3.8L) CS/4

LOT



For In Vitro Diagnostic Use

59°F (15°C)  86°F (30°C)

Date Received: _____ Date Opened: _____

10% Alcoholic Formalin is designed to rapidly fix tissues that are hard to fix, such as fatty tissues, breast and colon. It also functions well as a primary fixative, or on the tissue processor as a secondary fixative in the first station or after the initial fixative station(s).



DANGER - CONTAINS FORMALDEHYDE HAZARD STATEMENT(S):

Flammable liquid and vapor. Harmful if inhaled. Causes skin irritation. Causes serious eye damage. May cause an allergic skin reaction. May cause cancer. Causes damage to organs.

PRECAUTIONARY STATEMENTS - PREVENTION

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Use personal protective equipment as required. Use only outdoors or in a well-ventilated area. Wash face, hands and any exposed skin thoroughly after handling. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves. Do not breathe dust/fume/gas/mist/vapors. Do not eat, drink or smoke when using this product. Keep away from heat/sparks/ open flames/hot surfaces. — No smoking. Keep container tightly closed. Ground/ bond container and receiving equipment. Use explosion-proof electrical/ ventilating/lighting/equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

PRECAUTIONARY STATEMENTS - RESPONSE

IF EXPOSED: Call a POISON CENTER or doctor/physician.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. Immediately call a doctor.

IF ON SKIN: Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse. If skin irritation or rash occurs: Get medical advice/attention.

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or doctor/physician if you feel unwell.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL: Store locked up in a well-ventilated place at room temperature. Keep cool. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com **Made in USA**

MX1980 Alcoholic Formalin Gallons Case of 4 AW-2664753 Rev.B 2022-03-18