



10% Alcoholic Zinc Formalin

#FX1024 1 Gallon

#FX1025 5 Gallon Cube

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____

Buffered to a neutral pH, and formulated for rapid fixation of tissue, especially those with abundant fat present.



DANGER HAZARD STATEMENTS

Harmful if inhaled. Causes severe skin burns and eye damage. Causes serious eye damage. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction. May cause genetic defects. May cause cancer. May cause damage to organs (nose, lungs). May cause damage to organs (nose, lungs) through prolonged or repeated exposure. May damage fertility or the unborn child.

PRECAUTIONARY STATEMENTS - PREVENTION

Prevention: Avoid breathing dust, vapors. Use only outdoors or in a well-ventilated area. Do not breathe dusts or mists. Wash body thoroughly after handling. Wear protective gloves, protective clothing, eye protection and face protection. Wear eye protection, face protection. Wear NIOSH approved respiratory protection. Contaminated work clothing must not be allowed out of the workplace. Wear protective gloves. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Do not breathe dust, vapors. Do not eat, drink or smoke when using this product.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. **IF ON SKIN (or hair):** Take off immediately all contaminated clothing. Rinse skin with water/shower. Take off all contaminated clothing and wash it before reuse. Wash with plenty of water. If skin irritation or rash occurs, or if exposed or concerned: Get medical attention. **IF INHALED:** If breathing is difficult, remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: Call a doctor. **IF SWALLOWED:** Rinse mouth. Do NOT induce vomiting. **IN CASE OF FIRE:**

Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL: Store locked up. Dispose of contents/container to an approved waste disposal plant. **NOTICE:** Read Safety Data Sheet for additional information. **EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
cancerdiagnostics.com

Made in USA