



CoverSeal™-D

D-LIMONENE MOUNTING MEDIUM

#FXD178 500ML

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____

CoverSeal™-D is a unique non-toxic Toluene and Xylene-free permanent mounting medium that is compatible with all clearing agents and exhibits excellent consistency without application bubbles. CoverSeal™-D can replace any routine coverslipping medium. CoverSeal™-D is soluble in Xylene.



WARNING HAZARD STATEMENTS

Causes skin and serious eye irritation. May cause an allergic skin reaction. Flammable liquid and vapor.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Avoid breathing dust, vapors. Contaminated work clothing must not be allowed out of the workplace. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. If eye irritation persists: get medical attention. If exposed or concerned: get medical attention/call a doctor. Call a doctor if you feel unwell.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store in a well-ventilated place. Keep cool. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



CANCER
DIAGNOSTICS INC.

4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com **Made in USA**